Instruction

PERTUSSCAN PT IgG

ELISA for detection of IgG against *Bordetella pertussis* toxin (PT)
in serum

Microtitration strips (12x8) 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Doc. No. E-23-0207-02
October, 2013

PERTUSSCAN PT IgG

English:	page	2
Français:	page	17
Español:	página	32
Deutsch:	Seite	47
Italiano:	pagina	62
Dansk:	side.....	77
Norsk:	side.....	92
Svenska:	sida.....	107

PERT-PTG2

INTENDED USE OF THE TEST

Pertuscan PT IgG is an enzyme linked immunosorbent assay (ELISA) for the detection of IgG antibodies against Bordetella pertussis toxin (PT) in serum specimens. The test is intended for determination of current/recent infection.

The analysis should be performed by trained laboratory professionals.

“For in vitro diagnostic use”.

SUMMARY AND EXPLANATION

Pertussis, or whooping cough, is a contagious bacterial disease caused by the bacterium Bordetella pertussis.

Symptoms are initially mild (resembles those of a common cold), but will develop into severe coughing fits, which produce the namesake high-pitched "whoop" sound in infected babies and children when they inhale air after coughing. The coughing stage lasts for approximately six weeks before subsiding. Prevention via vaccination is of primary importance as treatment is of little clinical benefit to the person infected. Antibiotics, however, do decrease the duration of infectiousness and are thus recommended in many countries. The Pertussis Toxin (PT) is of significant importance for the pathogenesis of whooping cough. It is a real exotoxin responsible for many physiological, immunological and pharmacological effects. In addition the Pertussis toxin is highly specific for B. pertussis in contrast to some other exotoxins of the species Bordetella. Serology with determination of IgG against PT has proven especially useful for later diagnosis of older children and adults.

PRINCIPLE OF THE TEST

The assay consists of microtitre strips coated with highly purified PT as antigen. In infected individuals the antibodies produced against PT will bind to its antigen. Alkaline phosphatase conjugated anti-human IgG will bind to the human antibodies attached to the antigen. When the substrate is added a yellow colour will develop, the intensity of which (read in a photometer), is proportional to the amount of bound antibodies.

MATERIALS PROVIDED

One frame with microtitre strips (12x8) wells coated with PT antigen, packed in foil pouch with desiccant

Controls

PT IgG positive control. Ready-to-use

1 vial containing 1.5 mL, blue solution

PT IgG negative control (labelled Green). Ready-to-use

1 vial containing 4.0 mL, blue solution

Sample dilution buffer. Ready-to-use

2 vials containing 25 mL, blue solution

Wash buffer

1 vial containing 35 mL, 40 x concentrated

Conjugate

Anti-human IgG conjugate. Ready-to-use

2 vials containing 6.0 mL, orange solution

Substrate tablets

1 Strip containing 6 x 1 tablet

Substrate buffer

1 vial containing 30 mL

Stop solution

1 vial containing 30 mL

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.
- Incubation blocks or incubator (37°C).
- Non-metallic container for preparation of substrate solution.

REAGENTS

All reagents should have reached room temperature when used.

Sample dilution buffer. Contains phosphate buffered saline (PBS), 0.09% sodium azide. Ready-to-use.

Wash buffer, 40 x concentrated. 29% sodium chloride, 2% Tween 20 and 0.09% sodium azide. Dilute before use.

Controls

One PT IgG positive human serum and one PT IgG negative human serum. All sera are tested HBsAg and HIV- and HCV- antibody negative. Both controls are ready-to-use.

The conjugate contains goat-anti-human IgG conjugated with alkaline phosphatase. Ready-to-use.

Substrate tablets. P-Nitrophenyl phosphate, 5 mg/tablet.

Please note! Contact with metal should be avoided.

Substrate buffer. Contains 1M Diethanolamine and 0.5 mM Magnesium Chloride. Ready-to-use.

Stop solution. Contains 0.4 M sodium hydroxide, EDTA, carbonate buffer (pH >10). Do not allow the reagent to get into contact with the skin. Ready-to-use.

SAFETY AND PRECAUTIONS

General

Serum samples to be tested may contain infectious agents and should be handled accordingly.

Do not pipette by mouth, smoke, eat or drink in areas in which specimens or kit reagents are handled.

Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

ELISA plate

The microtitre strips should be handled with caution. To avoid contamination of wells never touch top of strips with fingers. Pipette carefully to avoid scratching the inside of the well. Check for air bubbles in wells after completed pipetting. If present, remove them by gently tapping the strip-holder.

Positive control

Contains human serum positive for PT IgG. Should be regarded as potentially infectious. Handle with caution. Tested negative for HIV- and HCV antibody and HbSAg.

Negative control

Contains human serum negative for PT IgG. Shall be regarded as potentially infectious. Handle with caution. Tested negative for HIV-and HCV-antibody and HbSAg.

Conjugate, sample diluent buffer


Contain biological material. Handle with caution.

Stop solution

Contains 0.4 M sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.

Warning

Contains: Sodium hydroxide

- H315: Causes skin irritation.
 H319: Causes serious eye irritation.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+352: IF ON SKIN: Wash with plenty of soap and water.
 P305+351 IF IN EYES: Rinse cautiously with water for several minutes.
 +338: Remove contact lenses, if present and easy to do. Continue rinsing.
 P332+313: If skin irritation occurs: Get medical advice/attention.
 P337+313 If eye irritation persists: Get medical advice/attention.

Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.

The Centers for Disease Control and Prevention and National Institutes of Health recommend that potentially infectious agents be handled at the Biosafety Level 2.

Substrate buffer

Contains 9,7% Diethanolamine

SUBSBUF

Warning

Contains: Diethanolamine

- H319: Causes serious eye irritation.
- P264: Wash hands thoroughly after handling.
- P280: Wear protective gloves/protective clothing/eye protection/face protection.
- P305+351 IF IN EYES: Rinse cautiously with water for several minutes.
- +338: Remove contact lenses, if present and easy to do. Continue rinsing.
- P337+313: If eye irritation persists: Get medical advice/attention.

STORAGE INSTRUCTIONS

General

Always store test kit and its contents at 2-8° C.
Test kit and its contents should have reached room temperature before use.

DO NOT USE KITS OR COMPONENTS THAT HAVE PASSED EXPIRY DATE.
DO NOT MIX REAGENTS WITH DIFFERENT BATCH NUMBERS.

Substrate solution

Prepared substrate solution should be used within the same day.
Keep in the dark at room temperature (20-23° C).

Wash solution

Use prepared wash solution within 3 weeks. Store at room temperature (20-23° C).

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Remove only the number of strips needed for testing, resealing the aluminium foil pouch carefully. Store at 2-8° C.

PREPARATION OF REAGENTS

All reagents should have reached room temperature when used.

Wash solution

Dilute 1 bottle of wash buffer (35 mL) to 1400 mL with distilled water. The prepared amount of wash buffer solution is sufficient for 12 microtitre strips and is stable for 3 weeks at 20-23° C.

Substrate solution

Prepare only the amount sufficient for one day use (100 µL/well).
1 substrate tablet dissolved in 5 mL substrate solution.

Please note! The substrate solution is stable for 1 day. Keep the prepared solution in the dark and not above room temperature.

WASHING INSTRUCTIONS

Always start the test procedure by washing each well four times.

Follow the washing instructions closely as incomplete washing may endanger the test result. Washing can be performed either by automatic microtitre plate washer or manually.

Washing with automatic microtitre plate washing equipment

When using automatic plate washing equipment, check that all wells can be aspirated completely and that the wash buffer is correctly dispensed reaching the rim of each well during each washing cycle. The washer should be programmed to execute 4 washing cycles. Continue immediately to the next reagent addition step.

Manual washing

1. Fill all the wells with 0.3 mL wash solution.
2. Let the filled wells stand for at least 30 seconds.
3. Empty the wells by turning the microtitre plate upside down followed by a firm short vertical movement.
4. This washing cycle (1, 2, and 3) should be carried out 4 times.
5. Place the inverted plate on absorbent paper towels and tap the plate firmly in order to remove residual wash solution in the wells.
6. Continue immediately to the next reagent addition step.

Do not scratch inside of well.

COLLECTION OF SAMPLES

Pertuscan PT IgG assay is for use with serum. Handle as if capable of transmitting infectious agents. Avoid using sera which are icteric, lipemic or hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The CLSI provides recommendations for storing blood specimens (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Paired samples

The most reliable method to diagnose a current/recent pertussis infection is by showing a significant titre rise (= a doubling of the absorbance value) in IgG between the acute and the convalescent sample.

An acute sample is defined as being drawn within:

2-3 weeks for unvaccinated children

<7 days for vaccinated children and adults, after onset of symptoms.

If the first sample has been drawn later than the indicated intervals, a significant rise in titre may not regularly be detected.

Single (convalescent) samples

Current/recent pertussis infection is also indicated by an elevated PT IgG (see cut-off under Interpretation) in convalescent samples, defined as a sample drawn after onset of symptoms

- ≥3 weeks for unvaccinated children
- ≥10 days for vaccinated children and adults

NB

Vaccination within 1 year prior to sampling renders the serological diagnosis uncertain.

TEST PROCEDURE

Always test each patient sample in duplicates.

Always use the supplied ready-to-use positive and negative controls in duplicates at each test run.

All reagents should be used at room temperature (20-23° C).

We recommend preparing a blank with only substrate solution. This is simply done by using an empty well to which is added substrate solution when adding substrate solution to the other wells.

1. Prepare wash solution. (See under "Preparation of Reagents").
2. Prepare the amount of substrate solution sufficient for one day use. (See under "Preparation of Reagents").
3. Dilute patient sera 1:500 (in duplicate) in test tubes in two steps.
Step I: 10 µL serum + 190 µL sample dilution buffer.
Step II: 15 µL serum from step I + 360 µL sample dilution buffer.
4. Place the appropriate number of strips into frame.
5. To maximize the activity of the PERTUSSCAN PT IgG assay, always start procedures by washing each well four times. (See under "Washing Instructions").
6. Pipette 100 µL of each ready-to-use control as follows.
PT IgG positive control.
Negative control (labelled Green).
DO NOT DILUTE THE CONTROLS.
7. Pipette 100 µL of each freshly prediluted patient serum in duplicate.
8. Incubate for 60 minutes at 37° C.
9. Wash each well four times. (See under "Washing Instructions").
10. Pipette 100 µL of conjugate into the wells.
DO NOT DILUTE THE CONJUGATE
11. Incubate for 60 minutes at 37° C.
12. Wash each well four times. (See under "Washing Instructions").
13. Pipette 100 µL substrate solution into all microtitre wells.
14. Incubate for 30 minutes at 37° C.
15. Pipette 100 µL stop solution into all microtitre wells.
16. Read absorbance values at 405 nm.

Note: Before calculation of test results, the OD value of the blank well should be subtracted from all other wells.

VALIDATION OF TEST

Calculate the mean of duplicate wells for each control.

	Target OD (at 405 nm)	Acceptable range
Positive Controls	1.00	0.77 – 1.43
Negative Control	-	<0.3

CALCULATION OF RESULTS

To correct for positive control that does not meet the target value 1.0, a calculation has to be made according to:

$$\frac{\text{OD}_{\text{Sample}}}{\text{OD}_{\text{Positive Control}}}$$

The corrected OD from the equation is used for interpretation of the result.

INTERPRETATION OF PERTUSSCAN PT IgG

(See also "Collection of Samples").

Calculate the mean for all duplicate patient samples.

Paired samples

A doubling of absorbance between first and second sample represents a significant rise of antibodies.

A significant rise of anti-PT-IgG indicates current/recent infection.

Single (convalescent) samples

A sample, drawn 4-6 weeks after onset of symptoms in an unvaccinated child and after ≥ 10 days of symptoms in a vaccinated child or an adult, showing an IgG titre:

≥ 1.5 strongly indicates a current/recent infection (specificity 99.1%) ≥ 1.0 indicates a current/recent infection (specificity 97.2%)

NB

- The absorbance value of one sample from one test run to another may slightly differ. This is mainly due to slight difference in test runs with regard to pipetting technique, temperature in incubator and laboratory, and time factor at incubations.
- Some specimens may not be repeatedly reactive. Common causes are:
 - a. improper washing
 - b. use of reagents which passed expiry date,
 - c. cross-contamination ("pick ups") from neighbour positive microtitre wells.
- A titre ≥ 1.5 is diagnostic for current/recent pertussis only if no pertussis vaccine has been given within 12 months prior to blood sampling.
- An absorbance value above 0.3 indicates presence of antibody
- Elderly individuals (>60 years) could have no elevated PT IgG response to a current/recent infection but may then be diagnosed by an elevated FHA IgA (≥ 0.6) response (Pertusscan TRO).
- If further confirmation of an elevated absorbance value of IgG to PT at ≥ 1.0 - < 1.5 is indicated then Pertusscan TRO can be used:
 - current pertussis will be confirmed if also IgA to PT and/or FHA is elevated
- The cut-off limit of 1.5 corresponds to ≥ 110 EU/mL against the US FDA standard Lot 3
The cut-off value of $\geq 1,0$ corresponds to ≥ 70 EU/mL
The cut-off value of 0.3 corresponds to 18 EU/mL

Pertusscan PT IgG test is controlled against NIBSC.

LIMITATIONS

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used as an adjunct to clinical symptoms and the results of other available tests.

PERFORMANCE CHARACTERISTICS

Table 1. **Clinical sensitivity and specificity.** A total of 159 frozen retrospective sera. The following table summarises the results:

Control and Disease groups	Total number	Positive Cut off OD ≥ 1.5	Positive Cut-off OD ≥ 1.0
PT IgG positive samples	52	43	51
Routine samples (non-respiratory infections)	107	1	3

Clinical sensitivity: 52 samples from patients diagnosed with Pertussis infection were tested for elevated PT IgG levels.

Cut-off OD ≥ 1.0 , sensitivity 51/52 = 98.1% 95% CI = 89.7% - 100%

Cut-off OD ≥ 1.5 , sensitivity 43/52 = 82.7% 95% CI = 69.7% - 91.8%

Clinical specificity: 107 Routine samples (non-respiratory infections) were tested.

Cut-off OD ≥ 1.0 , specificity 104/107 = 97.2% 95% CI = 92.0% - 99.4%

Cut-off OD ≥ 1.5 , specificity 106/107 = 99.1% 95% CI = 94.9% - 100%

The 95% confidence interval (CI) was calculated using the exact method.

Table 2. **Batch to batch variation** was determined by testing ten different samples in eight replicates on three different batches.

	1	2	3	4	5	6	7	8	9	10
Mean OD value	0.43	0.12	0.07	0.20	0.76	1	1.23	0.89	2.12	1.86
SD	0.02	0.01	0.01	0.01	0.07	0.03	0.12	0.06	0.04	0.13
% CV	3.8	8.6	9.2	6.9	8.7	2.6	9.9	6.6	2.1	6.9

Table 3. **Inter-assay precision** was determined by testing ten different samples in eight replicates at three separate occasions.

	1	2	3	4	5	6	7	8	9	10
Mean OD value	0.42	0.11	0.06	0.19	0.71	0.97	1.28	0.83	2.18	1.72
SD	0.01	0	0	0.01	0.01	0.02	0.02	0.03	0.05	0.02
% CV	1.5	2.7	4.6	3.6	1.7	1.6	1.7	3.8	2.4	0.9

Table 4. **Intra-assay precision** was determined by testing ten different samples in eight replicates at one occasion.

	1	2	3	4	5	6	7	8	9	10
Mean OD value	0.42	0.11	0.06	0.2	0.72	0.98	1.3	0.86	2.15	1.7
SD	0.03	0	0	0.01	0.04	0.05	0.07	0.03	0.1	0.1
% CV	6.1	4.0	5.3	4.9	4.8	5.8	4.8	3.9	3.5	5.8

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagent. Repeat test.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water solution used to prepare. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

REFERENCES

1. Cherry JD, et al. Defining Pertussis Epidemiology: Clinical, Microbiologic and Serologic Perspectives. *Ped Inf Dis J* 2005; 24:S25-34
2. Dalby T, et al. Pertussis serology: assessment of IgG anti-PT ELISA for replacement of the CHO cell assay. *APMIS* 2010; 118:968-972
3. Diez-Domingo J, et al. Incidence of pertussis in persons \leq 15 years of age in Valencia, Spain: seroprevalence of antibodies to pertussis toxin (PT) in children, adolescents and adults. *J Infection* 2004; 49:242-247
4. Granström G, Wretling B, Salenstedt C-R, Granström M. Evaluation of serologic assays for diagnosis of whooping cough. *J Clin Microbiol* 1988;1818-1823
5. Guiso N et al. What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* 2011;30:307–312
6. Hallander HO, et al. Seroprevalence of pertussis antitoxin (anti-PT) in Sweden before and 10 years after the introduction of a universal childhood pertussis vaccination program. *APMIS* 2009; 117:912-922
7. Menzies SL et al. Development and Analytical Validation of an Immunoassay for Quantifying Serum Anti-Pertussis Toxin Antibodies Resulting from *Bordetella pertussis* Infection. *Clin Vaccine Immunol* 2009;16:1781–1788
8. Riffelmann M, et al. Performance of Commercial Enzyme-Linked Immunosorbent Assays for Detection of Antibodies to *Bordetella pertussis* *J Clin Microbiol* 2010;48:4459-4463
9. Storsaeter J, Hallander HO, Gustafsson L, Olin P. Low levels of antipertussis antibodies plus lack of history of pertussis correlate with susceptibility after household exposure to *Bordetella pertussis*. *Vaccine* 2003; 21:3542-3549

Explanation of symbols.

	Batch code
	Catalogue number
	Use-by date
	Temperature limit
	Biological risks
	Consult instructions for use
	In vitro diagnostic medical device
	Manufacturer
	Warning
	Contains sufficient for 96 tests
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive.

Ag	Antigen (coated strips)
DIL	Sample dilution buffer
WASHBUF 40x	Wash buffer 40 x concentrated
CONTROL -	PT IgG negative serum
CONTROL + PT IgG	PT IgG positive serum
CONJ PT IgG	PT IgG conjugate
SUBSBUF	Substrate buffer
SUBS TAB pNPP	Substrate tablets
SOLN STOP	Stop solution (sodium hydroxide, 0.4 M)

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BIBLIOGRAPHIE

1. Cherry JD, et al. Defining Pertussis Epidemiology: Clinical, Microbiologic and Serologic Perspectives. *Ped Inf Dis J* 2005; 24:S25-34
2. Dalby T, et al. Pertussis serology: assessment of IgG anti-PT ELISA for replacement of the CHO cell assay. *APMIS* 2010; 118:968-972
3. Diez-Domingo J, et al. Incidence of pertussis in persons ≤ 15 years of age in Valencia, Spain: seroprevalence of antibodies to pertussis toxin (PT) in children, adolescents and adults. *J Infection* 2004; 49:242-247
4. Granström G, Wretling B, Salenstedt C-R, Granström M. Evaluation of serologic assays for diagnosis of whooping cough. *J Clin Microbiol* 1988;1818-1823
5. Guiso N et al. What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* 2011;30:307–312
6. Hallander HO, et al. Seroprevalence of pertussis antitoxin (anti-PT) in Sweden before and 10 years after the introduction of a universal childhood pertussis vaccination program. *APMIS* 2009; 117:912-922
7. Menzies SL et al. Development and Analytical Validation of an Immunoassay for Quantifying Serum Anti-Pertussis Toxin Antibodies Resulting from *Bordetella pertussis* Infection. *Clin Vaccine Immunol* 2009;16:1781–1788
8. Riffelmann M, et al. Performance of Commercial Enzyme-Linked Immunosorbent Assays for Detection of Antibodies to *Bordetella pertussis* *J Clin Microbiol* 2010;48:4459-4463
9. Storsaeter J, Hallander HO, Gustafsson L, Olin P. Low levels of antipertussis antibodies plus lack of history of pertussis correlate with susceptibility after household exposure to *Bordetella pertussis*. *Vaccine* 2003; 21:3542-3549

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REFERENSER

1. Cherry JD, et al. Defining Pertussis Epidemiology: Clinical, Microbiologic and Serologic Perspectives. *Ped Inf Dis J* 2005; 24:S25-34
2. Dalby T, et al. Pertussis serology: assessment of IgG anti-PT ELISA for replacement of the CHO cell assay. *APMIS* 2010; 118:968-972
3. Diez-Domingo J, et al. Incidence of pertussis in persons \leq 15 years of age in Valencia, Spain: seroprevalence of antibodies to pertussis toxin (PT) in children, adolescents and adults. *J Infection* 2004; 49:242-247
4. Granström G, Wretling B, Salenstedt CR, Granström M. Utvärdering av serologiska tester för diagnos av kikhosta. *J Clin Microbiol* 1988;1818-1823
5. Guiso N et al. What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* 2011; 30:307-312
6. Hallander HO, et al. Seroprevalence of pertussis antitoxin (anti-PT) in Sweden before and 10 years after the introduction of a universal childhood pertussis vaccination program. *APMIS* 2009; 117:912-922
7. Menzies SL et al. Development and Analytical Validation of an Immunoassay for Quantifying Serum Anti-Pertussis Toxin Antibodies Resulting from *Bordetella pertussis* Infection. *Glin Immunol Vaccine* 2009; 16:1781-1788
8. Riffelmann M, et al. Performance of Commercial Enzyme-Linked Immunosorbent Assays for Detection of Antibodies to *Bordetella pertussis* *J Clin Microbiol* 2010;48:4459-4463
9. Storsaeter J, Hallander HO, Gustafsson L, Olin P. Low levels of antipertussis antibodies plus lack of history of pertussis correlate with susceptibility after household exposure to *Bordetella pertussis*. *Vaccine* 2003; 21:3542-3549

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EURO DIAGNOSTICA AB