



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

Quick-DNA/RNA™ Pathogen Miniprep

Catalog Nos. R1042 & R1043

Highlights

- Spin-column purification of pathogen (virus, bacteria, protozoa) DNA/RNA from a wide variety of vectors (mosquitoes, fleas, ticks, *etc.*) and tissue types (mammals, birds, *etc.*) stored in DNA/RNA Shield™
- DNA/RNA is ready for any sensitive downstream applications (e.g., Next-Gen sequencing, RT/PCR, *etc.*)

Contents

Product Contents	1
Specifications	1
Product Description	2
Buffer Preparation	3
Protocol	3-4
Appendices	5
Ordering Information	6

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product please contact us.

Note: Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Product Contents

Quick-DNA/RNA™ Pathogen Miniprep (Kit Size)	R1042 (50 preps)	R1043 (200 preps)
DNA/RNA Shield™	50 ml	250 ml
Pathogen DNA/RNA Buffer¹	50 ml	100 ml
Proteinase K w/ Storage Buffer²	5 mg	20 mg
Zymo-Spin IIC Columns	50	200
Collection Tubes	50	200
Pathogen DNA/RNA Wash Buffer (concentrate)³	2x 6 ml	48 ml
DNase/RNase-Free Water	4 ml	2x 10 ml
Instruction Manual	1	1
ZR Bashing Bead™ Lysis Tubes (sold separately)	S6014-50 (1x 50 pack or 4x 50 pack)	

Storage Temperature – Store all kit components (i.e., buffers, columns) at room temperature.

¹ Add beta-mercaptoethanol to 0.5% (v/v) i.e., add 250 µl or 500 µl per 50 ml or 100 ml **Pathogen DNA/RNA Buffer**, respectively.

² Add 260 µl or 1,040 µl **Proteinase K Storage Buffer** per vial to reconstitute the lyophilized **Proteinase K**, 5 mg or 20 mg respectively. Vortex to dissolve and store frozen aliquots.

³ Add 24 ml of 100% ethanol (26 ml of 95% ethanol) to the 6 ml **Pathogen DNA/RNA Wash Buffer** concentrate (R1042) or 192 ml of 100% ethanol (204 ml of 95% ethanol) to the 48 ml **Pathogen DNA/RNA Wash Buffer** concentrate (R1043) before use.

Specifications

- **Sample Sources** – vectors (mosquitoes, fleas, ticks, other tough-to-lyse insects) and tissue types (animal tissue, plants, other hosts) processed and stored in DNA/RNA Shield™.
- **Sample Size** – ≤200 µl (or ≤10 mg).
- **Format** – 50 µg DNA/RNA binding capacity (spin-column), ≥25 µl elution volume
- **Size Range** – 50 nt to ~200 kb
- **Purity** – DNA/RNA is ready for any sensitive downstream applications (e.g., Next-Gen sequencing, RT/PCR, etc.)
- **Materials** (sold separately) –
ZR BashingBead Lysis Tubes; (S6014; 0.1/2.0 mm)

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

ZYMO RESEARCH CORP.


Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

Product Description


Quick-DNA/RNA™ Pathogen Miniprep kit is a spin-column based purification of pathogen (virus, bacteria, protozoa) DNA and RNA from a wide variety of vectors (mosquitoes, fleas, ticks, *etc.*) and tissue types (mammals, birds, *etc.*) collected, transported and stored in **DNA/RNA Shield™**.

The kit features a storage/lysis buffer system and can be combined with high density ZR BashingBead™ Lysis Tubes (optional) to facilitate complete homogenization of hard-to-lyse samples for efficient nucleic acid isolation. Small (>50 nt) and large (>200 kb) DNA and RNA are bound to **Zymo-Spin™ IIC Columns**, washed and then eluted.


The isolated high-quality nucleic acids are suitable for all downstream applications such as Next-Gen sequencing, hybridization-based and RT/PCR detection.



Inhibitor-free detection of West Nile Virus in mosquitoes



High-sensitivity detection of HIV-1 virus at low titer



(A) West Nile Virus (spiked-in mosquito homogenate),
(B) HIV-1 viral RNA particles (spiked-in plasma),
purified using the **Quick-DNA/RNA™ Pathogen** kit and
detected by RT-qPCR.

For **Assistance**, please
contact Zymo Research
Technical Support at
1-888-882-9682 or e-mail
tech@zymoresearch.com.

*optional

ZR BashingBead Lysis Tubes 2.0 mm + 0.1 mm
(catalog no. S6014)



ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

All steps should be performed at room temperature unless specified.

Buffer Preparation

- Add beta-mercaptoethanol 250 μ l or 500 μ l per 50 ml or 100 ml **Pathogen DNA/RNA Buffer**, respectively. (final concentration 0.5% (v/v)).
- Add 24 ml of 100% ethanol (26 ml of 95% ethanol) to the 6 ml **Pathogen DNA/RNA Wash Buffer** concentrate (R1042) or 192 ml of 100% ethanol (204 ml of 95% ethanol) to the 48 ml **Pathogen DNA/RNA Wash Buffer** concentrate (R1043) before use.
- Add 260 μ l or 1,040 μ l **Proteinase K Storage Buffer** per vial to reconstitute the lyophilized **Proteinase K**, 5 mg or 20 mg respectively. Vortex to dissolve and store frozen aliquots.

Protocol

- All steps should be performed at room temperature, unless specified.

Sample Preparation

1. Add **DNA/RNA Shield™** to the sample as recommended (table below) in a nuclease-free tube (not provided) and mix well:

		add DNA/RNA Shield™
		#R1100
Insects ¹ (tough-to-lyse, ticks, mosquitoes, fleas, deer flies, etc.)	≤ 10 mg	800 μ l
Tissue (mammals, birds, plants)		

For biological samples² collected and stored in **DNA/RNA Shield™** (i.e., #R1200, R1101, R1150, R1104, R1107, R1109, etc.), proceed to step 3 below.

2. Optional: Mechanical homogenization with the ZR BashingBead™ Lysis Tube³ and a high-speed cell disruptor⁴ is recommended for tough-to-lyse insects, animal tissue, plants, etc.

	homogenization time (high-speed)
insects (tough-to-lyse, ticks, mosquitoes, fleas, deer flies, etc.)	3-5 minutes
tissue (mammals, birds, plants)	30-60 seconds

3. To remove particulate debris or precipitation, centrifuge at 10,000-16,000 x g for 1 minute. Transfer up to 200 μ l of the cleared supernatant⁵ into a nuclease-free tube or well/plate (not provided).

Notes:

¹ See Appendices, page 5 for specific number of insect input.

² Plasma, serum, whole-blood, urine, fecal, swab, saliva, cell suspension, culture media, etc.

³ S6014-50, S6014-200

⁴ Required homogenization time will vary depending on the device and application. For high-speed cell disruptors (e.g., FastPrep® - 24, TerraLyzer™ Sample Processor or similar), samples can be processed in ≤ 5 minutes. For low-speed cell disruptors (e.g., Disruptor Genie™, or standard benchtop vortexes), processing can be ≤ 20 minutes long. See manufacturer's literature for operating information.

⁵ Up to 200 μ l liquid sample can be processed per prep.

DNA/RNA Purification

- Perform all steps at room temperature and centrifugation 10,000-16,000 x g for 30 seconds, unless specified.
4. Add 2 µl **Proteinase K** and mix well.
 5. Add 400 µl **Pathogen DNA/RNA Buffer**⁶ to each 200 µl sample⁵, mix well and incubate at room temperature for 5 minutes.
 6. Transfer the mixture into a **Zymo-Spin™ IIC Column**⁷ in a **Collection Tube** and centrifuge. Discard the flow-through.
 7. Add 500 µl **Pathogen DNA/RNA Wash Buffer**⁸ to the column and centrifuge. Discard the flow-through and repeat this step.
 8. Add 500 µl ethanol (95-100%) to the column and centrifuge for 1 minute to ensure removal of any residual ethanol. Discard the collection tube and carefully transfer the column into a new nuclease-free tube (not provided).
 9. Add 50 µl **DNase/RNase-Free Water** directly to the matrix of the column and centrifuge to elute the DNA/RNA.

Alternatively, for highly concentrated DNA/RNA use ≥25 µl elution volume.

The eluted DNA/RNA⁹ can be used immediately or stored frozen.

Notes:

⁵ Up to 200 µl liquid sample can be processed per prep.

⁶ To ensure efficient lysis and deproteinization, up to 5 volumes of Pathogen DNA/RNA Buffer can be used per 200 µl liquid sample.

⁷ Volume capacity of the column is 800 µl.

⁸ Before starting, add the appropriate amount of ethanol to the wash buffer concentrate, see Buffer Preparation, page 3.

⁹ It is recommended to titrate the DNA/RNA eluate for downstream applications (i.e., RT/PCR, etc.).

Appendices

Insect Samples: Recommended Input

Sample type	Maximum specimen input per 800 µl DNA/RNA Shield™
mosquito	≤ 50
tick	1 engorged of any species ≤ 5 flat adults ≤ 20 nymphs
flea	≤ 10
deer fly	1 adult

Ordering Information

Product Description	Catalog No.	Kit Size
Quick-DNA/RNA™ Pathogen Miniprep	R1042	50 Preps
	R1043	200 Preps
Quick-DNA/RNA™ Pathogen MagBead	R2145	96 Preps
	R2146	4x 96 Preps

For Individual Sale	Catalog No.	Amount
DNA/RNA Shield™	R1100-50	50 ml
	R1100-250	250 ml
Pathogen DNA/RNA Buffer	R1042-1-50	50 ml
	R1042-1-100	100 ml
Proteinase K w/ Storage Buffer Set	D3001-2-5	5 mg
	D3001-2-20	20 mg
Zymo-Spin™ IIC Columns	C1011-50	50 pack
	C1011-250	200 pack
Collection Tubes (2 ml)	C1001-50	50 pack
	C1001-500	500 pack
Pathogen DNA/RNA Wash Buffer (concentrate)	R1042-2-6	12 ml
	R1042-2-48	48 ml
DNase/RNase-Free Water	W1001-30	30 ml
	W1001-100	100 ml
	W1001-200	200 ml
ZR BashingBead™ Lysis Tubes (0.1 mm & 2.0 mm)	S6014-50	50 pack

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

DNA & RNA MADE SIMPLE



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com