

ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

EZ RNA MethylationTM Kit

Catalog Nos. **R5001 & R5002**

Highlights

- Fast and reliable bisulfite conversion of RNA for methylation analysis.
- Specifically optimized for complete conversion of non-methylated cytosine in RNA.
- Ideal for all RNA inputs.

Contents

| | |
|--------------------------------------|------|
| Product Contents | 1 |
| Product Specifications..... | 1 |
| Introduction to RNA Methylation..... | 2 |
| Product Description..... | 3 |
| Reagent Preparation | 4 |
| Protocol..... | 4 |
| Appendices | 5, 6 |
| Frequently Asked Questions | 7 |
| Ordering Information | 8 |
| List of Related Products | 9 |

For Research Use Only

Ver. 1.0.4

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-949-266-9452 • Web: www.zymoresearch.com • E-mail: info@zymoresearch.com

Product Contents:

| | R5001 | R5002 | Storage Temperature |
|--|--------------|--------------|----------------------------|
| EZ RNA Methylation™ Kit | 50 rxns. | 200 rxns. | |
| RNA Conversion Reagent | 5 tubes | 20 tubes | Room Temp. |
| RNA Binding Buffer | 25 ml | 100 ml | Room Temp. |
| RNA Wash Buffer¹ (concentrate) | 12 ml | 48 ml | Room Temp. |
| RNA Desulphonation Buffer | 10 ml | 40 ml | Room Temp. |
| DNase/RNase-Free Water | 1 ml | 4 ml | Room Temp. |
| Zymo-Spin™ IC Columns | 50 columns | 200 columns | Room Temp. |
| Collection Tubes | 50 tubes | 200 tubes | Room Temp. |
| Instruction Manual | 1 | 1 | – |

Note - Integrity of kit components is guaranteed for one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

¹ Add 48 ml 100% ethanol (52 ml of 95% ethanol) to the 12 ml **RNA Wash Buffer** concentrate (R5001) or 192 ml 100% ethanol (208 ml of 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate (R5002) before use.

Specifications:

- **RNA Input:** Samples containing 32 ng - 3 µg of DNA-free RNA. For optimal results, the amount of input RNA should be between 0.5 - 1 µg.
- **Conversion Efficiency:** > 99% of non-methylated C residues are converted to U with > 99% protection of 5-methylcytosine.
- **RNA Recovery:** > 80%

Note: For purification of high-quality DNA-free RNA, we recommend the **Quick-RNA™** (R1050) or **Direct-zol™** (R2050) purification kits, see page 9 for details. DNase I treatment of RNA samples is recommended.

Note: m4mC is also partially resistant to conversion with bisulfite, however, in comparison with 5-mC, m4mC may be more easily converted to U during the procedure.


Use of Methylation Specific PCR (MSP) is protected by US Patents 5,786,146 & 6,017,704 & 6,200,756 & 6,265,171 and International Patent WO 97/46705. No license under these patents to use the MSP process is conveyed expressly or by implication to the purchaser by the purchase of this product.

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

Introduction to RNA Methylation:

Although the majority of nucleic acid modification research involves 5-methylcytosine in DNA, RNA is also extensively modified. In fact, there exist over a hundred modifications to RNA (1). 5-methylcytosine (5-mC) is present in RNA, and methylation is a common and naturally-occurring event in the RNA from both prokaryotic and eukaryotic organisms (2, 3). However, the function of RNA methylation remains unknown. Some reports describe a role for RNA methylation in translational regulation (4, 5), while others support a role for methylation in regulating RNA stability (6, 7) or the facilitation of RNA structure formation (8, 9).

The ability to detect and quantify 5-methylcytosine in RNA efficiently and accurately has been troublesome due to the inability of RNA to withstand the pH and temperatures used in the standard workflow for bisulfite conversion of DNA. Zymo Research offers a solution to these problems with the **EZ RNA Methylation™ Kit** that has been optimized and validated for bisulfite conversion of RNA. This technique involves treating RNA with a unique bisulfite conversion reagent, which converts non-methylated cytosines into uracil while preserving the integrity of the RNA. Methylated cytosines remain unchanged during the treatment. After performing bisulfite treatment, the methylation profile of the RNA can be determined using techniques like RT-PCR followed by DNA sequencing (see figure on page 3).



The Chemistry of Bisulfite Conversion

References:

1. Cantara WA, *et al.* Nucleic Acids Res. 2011; 39: D195-201.
2. Motorin Y, Helm M. Nucleic Acids Res. 2010; 38(5): 1415-1430.
3. Squires JE, Preiss T. Epigenomics. 2010; 2(5): 709-715.
4. Chow CS, *et al.* ACS Chem Biol. 2007; 2(9): 610-619.
5. Baudin-Baillieu A, *et al.* Nucleic Acids Res. 2009; 37(22): 7665-7677.
6. Alexandrov A, *et al.* Mol Cell. 2006; 21(1): 87-96.
7. Schaefer M, *et al.* Genes Dev. 2010; 24(15): 1590-1595.
8. Helm M. Nucleic Acids Res. 2006; 34(2): 721-733.
9. Motorin Y, Helm M. Biochemistry. 2010; 49(24): 4934-4944.


Selected EZ RNA Methylation™ Kit Citation:

Amort T, *et al.* RNA Biol. 2013; 10(6): PMID 23595112.

For bisulfite treatment of DNA for methylation analysis, see the comprehensive line of **EZ DNA Methylation™** products from Zymo Research. (D5001-D5047)

Product Description:

The **EZ RNA Methylation™ Kit** features rapid and reliable bisulfite treatment and conversion of cytosines in RNA for methylation analysis. The kit streamlines the three-step process for complete conversion of cytosine into uracil. RNA denaturation and bisulfite conversion processes are combined into a single step. No buffer preparation is necessary. The **RNA Conversion Reagent** is provided ready-to-use: simply add the reagent to an RNA sample and incubate as indicated. Also, innovative in-column desulphonation technology eliminates messy precipitation steps, ensuring researchers obtain consistent results. The product has been designed to minimize template degradation, loss of RNA during treatment and clean-up, and to provide complete conversion of cytosine for accurate methylation analysis. Recovered RNA is ideal for RT-PCR, sequencing, library preparation and Next-Gen sequencing.



Original RNA with Methylated CpG:


AUCCUUC^mGAUGUCGG

Converted RNA:

AUUUUUC^mGAUGUUUGG

cDNA Sequence:

ATTTTTC^mGATGTTGG



Sequencing results following bisulfite treatment. RNA with methylated C (5-mC) at nucleotide position #7 was processed using the **EZ RNA Methylation™ Kit**. The recovered RNA was amplified by RT-PCR and then cloned and sequenced. The methylated cytosine at position #7 remained intact while the non-methylated cytosines at positions #3, 4, and 13 were completely converted into uracil (post-bisulfite treatment) and detected as thymine following RT-PCR and sequencing.

Reagent Preparation:

- **Preparation of RNA Wash Buffer:** Add 48 ml of 100% ethanol (52 ml of 95% ethanol) to the 12 ml **RNA Wash Buffer** concentrate or 192 ml of 100% ethanol (208 ml of 95% ethanol) to the 48 ml **RNA-Wash Buffer** concentrate before use.

Protocol for Bisulfite Conversion of RNA:

1. Add 130 μ l of **RNA Conversion Reagent** to 20 μ l of RNA sample in a PCR tube. Mix the sample by flicking the tube or pipetting up and down, then centrifuge briefly to ensure there are no droplets in the cap or sides of the tube.

Note: If the sample volume is less than 20 μ l, compensate with DNase/RNase-Free Water.

2. Place the PCR tube(s) in a thermal cycler and perform the following steps:

1. 70°C for 5 minutes
2. 54°C for 45 minutes
3. 4°C up to 20 hours

Note: The 4°C storage step is *optional*.

3. Place a **Zymo-Spin™ IC Column** into a **Collection Tube** and add 250 μ l of **RNA Binding Buffer** to the column.
4. Load the sample (from Step 2) into the **Zymo-Spin™ IC Column** containing the **RNA Binding Buffer** and mix by pipetting up and down.
5. Add 400 μ l of 95-100% ethanol to the sample-**RNA Binding Buffer** mixture in the column. Close the cap and immediately mix by inverting the column several times.
6. Centrifuge at full speed ($\geq 10,000 \times g$) for 30 seconds. Discard the flow-through.
7. Add 200 μ l **RNA Wash Buffer** to the column and centrifuge at full speed for 30 seconds.
8. Add 200 μ l of **RNA Desulphonation Buffer** to the column and let stand at room temperature (20°C – 30°C) for 30 minutes. After the incubation, centrifuge at full speed for 30 seconds. Discard the flow-through.
9. Add 400 μ l **RNA Wash Buffer** to the column and centrifuge at full speed for 30 seconds. Repeat the wash step with an additional 400 μ l **RNA Wash Buffer**. Discard the flow-through.
10. Centrifuge the **Zymo-Spin™ IC Column** in the emptied **Collection Tube** at full speed for 2 minutes. Remove the **Zymo-Spin™ IC Column** carefully from the **Collection Tube** and transfer it into an RNase-free Tube.
11. Add $\geq 10 \mu$ l of **DNase/RNase-Free Water** directly to the column matrix and let stand for 1 minute at room temperature. Centrifuge at full speed for 30 seconds. The eluted RNA can be used immediately or stored at -20°C for up to 3 months. For long-term storage, keep at or below -70°C.

Note: The elution volume can be $> 10 \mu$ l depending on the requirements of your experiments.

Appendix I: Using 28S rRNA as a Positive Control

We recommend using 28S ribosomal RNA (*H. sapiens*) as a positive control for RNA methylation analysis, as the C at position 4447 (GenBank accession # NR_003287) is generally 100% methylated. Total RNA from cells or tissues (i.e. HeLa, HCT116, keratinocytes, brain tissue, liver tissue, etc.) can be used directly for the bisulfite conversion. The following sequence is the 28S rRNA region amplified (post-conversion) using the primer set indicated below.

Original (non-converted):

```
4321 -----ggg gccucacgau ccuucugacc uuuuggguuu uaagcaggag gugucagaaa aguuaccaca
4391 gggauaacug gcuuguggcg gccaacgguu cauagcgacg ucgcuuuuug auccuuCgau gucggcucuu
4461 ccuaucuuug ugaagcagaa uucaccaagc guuggauugu ucaccacua auagggaaagc ugagcugg--
```

Bisulfite-Converted:

```
4321 -----ggg guuuuugau uuuuuugau uuuuggguuu uaaguaggag guuuuagaaa aguuuuuuaa
4391 gggauuuug guuuguggug guuuaguuu uauagugaug uuguuuuuug auuuuuCgau guugguuuuu
4461 uuuuuuuug ugaaguagaa uuuuuuuagu guuggauugu uuuuuuuuaa auagggaaug ugaguugg--
```

H 28SF primer: 5' -GGGGTTTTAYGATTTTTTTGATTTTTTTGGG-3'

H 28SR primer: 5' -CCAACACR'TTCCCTATTAATAAATAAAC-3'

Representative sequencing data obtained using 28S rRNA: The RT-PCR sequencing results of 10 clones (below) were obtained using bisulfite-converted total RNA extracted from HeLa cells. Underlined **C** represents 5-mC, highlighted **C** represents non-converted cytosine, *italics* are primer regions. The original, non-converted RNA sequence with non-methylated **C** highlighted is shown below the converted cDNA sequencing results for comparison.

Conversion Efficiency (C to T): C: 99.5%

```

1
HeLa01 - GGGGTTTTATGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa02 - GGGGTTTTATGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa03 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa04 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa05 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa06 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa07 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa08 - GGGGTTTTACGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa09 - GGGGTTTTATGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA
HeLa10 - GGGGTTTTATGATTTTTTGATTTTTGGGTTTTAAGTAGGAGGTGTAGAAAAGTTATTATAGGGATAATTGGTTTGTGGTGGTTAAGTGTTTATAGTGA

Orig. - GGGGCCUCA1CGAUCCUUCUGACC2UUUUGGUUUUAAGCAGGAGGUGTCAGAAAAGUUA3CCAGGGAAUACUGGC4UUGUGCGGCCAAG5GUUCAUAG6CA

102
HeLa01 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa02 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa03 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa04 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa05 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa06 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa07 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa08 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa09 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG
HeLa10 - TGTGTTTTTTGATTTTTCGATGTTGGTTTTTTTTATTATTGTGAAGTAGAATTTATTAAGTGTGGATTGTTTTATTTATTAATAGGGAATGTGAGTTGG

Orig. - CGUCC1UUUUUGAUCC2UUCGAUGUCCG3GUUUCCUAUCAUUGUGAAGCAGAAUUCACCAAG4GUUGGAUUGUUC5CCCAUAAUAGGGAA6GUGAGCUGG

```


Frequently Asked Questions:

Q: Should the input RNA be dissolved in TE, water, or some other buffer prior to its conversion?

A: Water is recommended.

Q: Should the RNA be DNase I treated?

A: Yes, it is recommended to treat RNA samples with DNase I (e.g., E1009) prior to the bisulfite conversion.

(For clean-up of DNase I-treated RNA use the RNA Clean & Concentrator™, R1015)

Q: At what temperature and for how long can converted RNA be stored?

A: The sample should be stored at $\leq -20^{\circ}\text{C}$ whenever possible and freeze-thaw cycles should be minimized. The quality of the RNA should remain relatively unchanged for up to 3 months. For long term storage samples should be kept at $\leq -70^{\circ}\text{C}$.

Q: Which Taq polymerase(s) do you recommend for PCR amplification of cDNA generated from bisulfite-converted RNA?

A: We recommend a “hot start” DNA polymerase (e.g., ZymoTaq™, E2001).

Q: What RNA purification methods do you recommend?

A: For RNA purification from cells or soft tissues use Quick-RNA™ kits (R1054). For samples in Tri-Reagent® or similar, the Direct-zol™ kits (D2050, D2051) are recommended. Both technologies allow for total RNA recovery (including small RNAs) and facilitate on-column DNase I treatment.

Ordering Information:

| Product Description | Catalog No. | Kit Size |
|-------------------------|-------------|-----------|
| EZ RNA Methylation™ Kit | R5001 | 50 rxns. |
| EZ RNA Methylation™ Kit | R5002 | 200 rxns. |

| For Individual Sale | Catalog No. | Amount(s) |
|--------------------------------|-------------|-------------|
| RNA Conversion Reagent | R5001-1-1 | 1 tube |
| RNA Binding Buffer | R1013-2-25 | 25 ml |
| | R1013-2-100 | 100 ml |
| RNA Wash Buffer (concentrate) | R1003-3-6 | 6 ml |
| | R1003-3-12 | 12 ml |
| | R1003-3-24 | 24 ml |
| | R1003-3-48 | 48 ml |
| RNA Desulphonation Buffer | R5001-3-10 | 10 ml |
| | R5001-3-40 | 40 ml |
| DNase/RNase-Free Water | W1001-1 | 1 ml |
| | W1001-4 | 4 ml |
| | W1001-6 | 6 ml |
| | W1001-10 | 10 ml |
| Zymo-Spin™ IC Columns (capped) | C1004-50 | 50 columns |
| | C1004-250 | 250 columns |
| Collection Tubes | C1001-50 | 50 tubes |
| | C1001-500 | 500 tubes |
| | C1001-1000 | 1,000 tubes |

Popular RNA Products from Zymo Research

THE
Epigenetics
COMPANY™

| Product | Description | Prep/Format | Catalog |
|--|---|--|-------------------------|
| RNA Clean-Up | | | |
| RNA Clean & Concentrator™-5 | Cleanup and concentration of modified, labeled, impure, diluted, DNase treated RNA (≥ 17nt) and purification of RNA from aqueous phase of organic extracts. (DNase I included with R1013 and R1014.) | 50/column 200/column | R1015 R1016 |
| RNA Clean & Concentrator™-25 | | 50/column 100/column | R1017 R1018 |
| RNA Clean & Concentrator™-100 | | 25/column | R1019 |
| ZR-96 RNA Clean & Concentrator™ | | 2x 96/plate | R1080 |
| DNA-Free RNA Kit™ | | 50/column 200/column | R1013 R1014 |
| Oligo Clean & Concentrator™ | Cleanup and concentration of RNA and/or DNA oligos. Good for clean-up of miRNAs and siRNAs. | 50/column 200/column | D4060 D4061 |
| ZR-96 Oligo Clean & Concentrator™ | | 2x 96/plate 4x 96/plate | D4062 D4063 |
| ssDNA/RNA Clean & Concentrator™ | | 20/column 50/column | D7010 D7011 |
| Zymoclean™ Gel RNA Recovery Kit | Recovery of RNA from agarose gels. | 50/column | R1011 |
| ZR small-RNA™ PAGE Recovery Kit | Small RNA (> 17nt) from polyacrylamide gels. | 20/column | R1070 |
| OneStep™ PCR Inhibitor Removal Kit | Removal of polyphenolics, humic/fulvic acids, tannins, melanin etc. from RNA. | 50/column | D6030 |
| OneStep™-96 PCR Inhibitor Removal Kit | | 2x 96/plate | D6035 |
| RNA from Samples in Trizol®, TRI Reagent®, etc. | | | |
| Direct-zol™ RNA MiniPrep | RNA (>17 nt) from TRI Reagent®, TRIzol®, and all other acid-guanidinium-phenol based reagents without phase separation. DNase I included. | 50/column 200/column | R2050 R2052 |
| Direct-zol™ RNA MiniPrep w/ TRI Reagent® | | 50/column w/ (50 ml) 200/column w/ (100 ml) | R2051 R2053 |
| Direct-zol™-96 RNA | | 2x 96/plate 4x 96/plate | R2054 R2056 |
| Direct-zol™-96 RNA w/ TRI Reagent® | | 2x 96/plate w/ (200 ml) 4x 96/plate w/ (400 ml) | R2055 R2057 |
| Direct-zol™-96 MagBead RNA | | 2x 96/plate 4x 96/plate 8x 96/plate | R2100 R2102 R2104 |
| Direct-zol™-96 MagBead RNA w/ TRI Reagent® | 2x 96/plate w/ (200 ml) 4x 96/plate w/ (400 ml) 8x 96/plate w/ (800 ml) | R2101 R2103 R2105 | |
| RNA from Cells and Tissue | | | |
| Quick-RNA™ MicroPrep | Total RNA from cells and tissue. DNase I included. | 50/column 200/column | R1050 R1051 |
| Quick-RNA™ MiniPrep | | 50/column 200/column | R1054 R1055 |
| Quick-RNA™ MidiPrep | | 25/column | R1056 |
| ZR-96 Quick-RNA™ | | 2x 96/plate 4x 96/plate | R1052 R1053 |
| ZR-Duet™ DNA/RNA MiniPrep | Parallel purification of DNA/RNA from cells. | 50/column | D7001 |
| Pinpoint™ Slide RNA Isolation System Kit I | RNA from fresh/frozen tissue sections. | 50/column | R1003 |
| Pinpoint™ Slide RNA Isolation System Kit II | RNA from paraffin-embedded (FFPE) tissue. | 50/column | R1007 |
| RNA from Biological Liquids | | | |
| ZR Viral RNA Kit™ | RNA (DNA) from body fluids (plasma, serum, CSF, urine). | 50/column 200/column | R1034 R1035 |
| ZR-96 Viral RNA Kit™ | | 2x 96/plate 4x 96/plate | R1040 R1041 |
| ZR Viral DNA/RNA Kit™ | | 25/column 100/column | D7020 D7021 |
| ZR Whole-Blood RNA MiniPrep™ | RNA from whole blood or partitioned blood. | 50/column 100/column | R1020 R1021 |
| ZR Urine RNA Isolation Kit™ | Cellular and endosomal RNA from urine. | 20/column 50/column | R1038 R1039 |
| RNA from Tough-to-Lyse Samples | | | |
| ZR Fungal/Bacterial RNA MicroPrep™ | RNA from bacteria, yeast, fungi; BashingBead™ lysis. | 50/column | R2010 |
| ZR Fungal/Bacterial RNA MiniPrep™ | | 50/column | R2014 |
| ZR Plant RNA MiniPrep™ | RNA from leaves, stems, buds, flowers, fruits, seeds, etc; BashingBead™ lysis, RT/PCR inhibitor removal. | 50/column | R2024 |
| ZR Tissue & Insect RNA MicroPrep™ | RNA from insect, arthropod specimen and small tissue samples; BashingBead™ lysis. | 50/column | R2030 |
| ZR Soil/Fecal RNA MicroPrep™ | RNA from soil, sludge, sediment, feces. | 50/column | R2040 |
| YeaStar RNA Kit™ | RNA from yeast strains susceptible to Zymolyase. | 50/column | R1002 |
| DNA and RNA Sample Preservation and Storage | | | |
| DNA/RNA Shield™ | Cells, biological liquid, tissue storage and DNA/RNA purification. | 50 ml 250 ml | R1100-50 R1100-250 |
| DNA/RNA Shield™ w/ Quick-RNA™ MiniPrep | | 50/column w/ (50 ml) | R1100 |
| Enzymes and Markers | | | |
| DNase I w/ 10X Reaction Buffer | Lyophilized | 250 U | E1009 |
| ZR small-RNA™ Ladder | ssRNA (17, 21, 25, 29 nt) | 10 µg | R1090 |

ZYMO RESEARCH CORP.

Toll Free: 1-888-882-9682 • Fax: 1-949-266-9452 • Web: www.zymoresearch.com • E-mail: info@zymoresearch.com