



Antibodies for Simple Western™

Certified Antibodies for ProteinSimple® Instruments:






- Guaranteed Simple Western performance
- · Largest selection of certified antibodies in the market with over 1,000 and growing
- Optimized antibody and lysate concentrations to save you time
- Antibodies receive a "Badge of Honor" which is easily visible in online search
- Test an antibody not yet certified in Simple Western and become eligible for rewards


See data examples below:

Traditional Western




Detection of Human Cytokeratin 19 by Western blot. Western blot shows detection of Cytokeratin 19 in untreated MCF-7 Lysate. The nitrocellulose membrane was probed with a Rabbit Anti-Human/Cytokeratin 19 Polyclonal Antibody (Catalog # NBP1-78278, Novus Biologicals). A specific band was detected for Cytokeratin 19 at approximately 40 kDa (as indicated).

Simple Western




Detection of Human Cytokeratin 19 by Simple Western. Simple Western lane view shows a specific band between 40-45 kDa for Cytokeratin 19 tested in 0.05 mg/mL of untreated MCF-7 Lysate using a Rabbit Anti-Human/Cytokeratin 19 Antibody (Catalog # NBP1-78278, Novus Biologicals). This experiment was performed under reducing conditions using the 12-230 kDa separation system.

Traditional Western



Detection of Human/Mouse/Rat HSP60 by Western blot. Western blot shows lysates of Jurkat human acute T cell leukemia cell line, MCF 7 human breast cancer cell line. NIH 3T3 mouse embryonic fibroblast cell line, and Nb2 11 rat lymphoma cell line. PVDF membrane was probed with 0.1 µg/mL of Human/Mouse/Rat HSP60 Antigen Affinity-purified Polyclonal Antibody (Catalog # AF1800) followed by HRP-conjugated Anti-Rabbit IgG Secondary Antibody (Catalog # HAF008, R&D Systems). A specific band was detected for HSP60 at approximately 62 kDa (as indicated).

Simple Western



Detection of Human and Rat HSP60 by Simple Western. Simple Western lane view shows lysates of Jurkat human acute T cell leukemia cell line and NIH 3T3 mouse embryonic fibroblast cell line, loaded at 0.2 mg/mL. A specific band was detected for HSP60 at approximately 60 kDa (as indicated) using 0.5 µg/mL of Rabbit Anti-Human/Mouse/Rat HSP60 Antigen Affinity-purified Polyclonal Antibody (Catalog # AF1800, R&D Systems). This experiment was conducted under reducing conditions and using the 12-230 kDa separation system.

Popular Antibodies

Snapshot of certified antibodies from both Novus Biologicals and R&D Systems

HDAC6
HGF R/c-MET
HIF-1 alpha
HIF-2 alpha
His Tag
Histone H4
HSP27
HSP60
HSP70/HSPA1A
IKK alpha
IKK beta
LC3
LDL Receptor
LOX
MMR/CD206
MTH1
NBR1
NBS1
Nestin
NeuN
NuMA
OXPAT
p38
p53
p62/SQSTM1
p70 S6 Kinase
p73
PCNA
PDZK1
Perilipin 2

phospho-p53
phospho-STAT1
phospho-STAT6
PINK1
PRAME
PTK7/CCK4
RBFOX3/NeuN
S100A7/Psoriasin
S100A9
SCP1
Sirtuin 1/SIRT1
SNAIL
SNX27
S0X2
SOX9
SR-BI
SREBP1
STAT5b
Survivin
TERT
Tie-2
TMEM26
TRIF
Tyrosine Hydroxylase
UCP1
Vimentin
VPS34
WT1
хСТ
ZEB1



LEARN MORE