


Protocol for IRODORI (Multi-color Western blot system)

Example of IRODORI



In the example 4 antigens (vimentin, GST, tubulin, actin) were detected using each antigen-specific antibodies which is attached to fluorescent MAD reagent labeled with 4 different Cy dyes.

Materials

fMAD (fluorescent MAD reagent with Cy dye-labeled):

4 different MAD types are available (Cy 2, Cy3, Cy4, Cy7).

BNC-ZZ (non labeled)

5 mM BS3 dissolved in PBS

10 mM glycine

Skim milk

TBST (0.02M Tris-HCl, 0.14M NaCl, 0.05% Tween 20, pH 7.4)

Protein A Sepharose

PBS (pH 7.4)

Antibodies

1. Preparation of antibody-fMAD complex (fMAD complex)

For multicolor western, prepare multiple complexes made of different antibody and labeled with Cy dyes of different fluorescence properties.

- 1) Dissolve 5 μ g of fMAD in adequate amount of PBS. To the fMAD solution, add 5 μ L of 1mg/mL of Antibody solution*, and 5 μ L of 5mM BS3 solution. Adjust the

volume to be 50 μL .

*: If the concentration of antibody solution less than 1mg/ml, increase the volume of antibody solution so that the amount of antibody should be 5 μg .

- 2) Incubate the solution for 30min at RT, and the at the end of incubation add 0.5 μL of 10mM glycine to stop the reaction
- 3) Add all the solution to 100 μL of Protein A Sepharose resin, incubate for 15min at RT, and centrifuge at 10000 rpm for 2 min at 4 .
- 4) Transfer 50 μL of the supernatant add 5 μg of BNC-ZZ , and incubation for 30 min at RT.
- 5) The antibody-fMAD complex is now ready.

2. Multi-color western

- 1) Prepare protein-transferred membrane by standard method.
- 2) Block the membrane by incubating with Blocking buffer (0.05% Skim milk dissolved in TBST) for 1hr.
- 3) When using multiple antibody-fMAD complexes, combine them. Then, dilute the antibody-fMAD complex solution to 500 μL by Blocking buffer .
- 4) Incubate the membrane with the diluted antibody-fMAD complex solution for 1hr at RT.
- 5) Take the fluorescent (pseudo-color) image of the membrane using corresponding excitation and emission wave lengths of each antibody-fMAD complex.
- 6) Combine the captured pseudo-color image.

Reference

Fluorophore-labeled nanocapsules displaying IgG Fc-binding domains for the simultaneous detection of multiple antigens

M.Iijima, T. Matsuzaki, N. Yoshimoto, T. Niimi, K. Tanizawa, S.Kuroda (2011)
Biomaterials 32, 9011-9020