


Quick-Start Protocol

Revision 1.0 Aug 2023

This protocol is designed to extract genomic DNA from fecal samples.

Notes before starting:

- Buffer SF2 needs to be stored at 2-8 °C upon reception.
- This protocol requires the use of a centrifuge capable of generating at least 15,000 g.
- For faster processing, pre-position the plasticwares used during the extraction as depicted below.



Column and sample preparation

Remove inhibitor

Bind

Wash

Elute

1. Add 200 µL of Equilibration Buffer to Column SF membranes to ensure its performance. Wait at least 1 min and centrifuge for 10 sec @ maximum speed. Transfer the Column SF into a new Collection tube (provided).

2. Weigh up to 250 mg of the feces and add it to a Lysing Matrix YB tube.

3. Add 900 µL of Buffer SF1. Homogenize using Fastprep 5 m/s for 35 sec or vortex at 2500-3000 rpm for 20 min, centrifuge for 2 min @ maximum speed.

4. During the centrifugation, Add 200 µL of Buffer SF2 into a new 2 mL Tube (provided). Transfer the supernatant (~500-700 µL) while avoiding the pellet, vortex for 1 sec and centrifuge for 2 min @ maximum speed.

5. Transfer the supernatant (~600-800 µL) into a 2 mL Tube (provided). Add 1 volume of Buffer SF3, vortex for 1 sec.

Note: if the supernatant and Buffer SF3 lysate is cloudy or include debris, centrifuge for 1 min @ maximum speed prior to binding to Column SF and transfer the supernatant as described below.

6. Apply ~750 µL of the lysate to Column SF, centrifuge for 10 sec @ 15,000 g and discard the flow-through. Repeat the process until all the lysate has passed through.

7. **1st wash.** Transfer Column SF into a new Collection tube (provided). Add 700 µL of Buffer SF4 to the center of the column, centrifuge for 10 sec @ 15,000 g. Discard the flow-through and place Column SF back into the same Collection tube.

8. **2nd wash.** Add 700 µL of Buffer SF5 to the center of the column and centrifuge for 30 sec @ ≥15,000 g.

9. **Column drying.** Transfer Column SF into a new Collection tube (provided), centrifuge for 1 min @ maximum speed.

10. **Elution.** Transfer the Column SF into a new Elution tube (provided). Add 100 µL of Buffer SF6 directly to the column membrane, wait at least 1 min and centrifuge for 1 min @ ≥15,000 g. The DNA sample is now ready for downstream applications.

Note: For more concentrated sample, elute with 50 µL of Buffer SF6.

SPINeasy™ DNA Pro Kit for Feces


Cat. No.: 116547050 (50 PREPS)/116547000 (5 PREPS)



Revision 1.0 Aug 2023

• Quick-Start Protocol

Column preparation
Sample preparation
Remove inhibitor
Bind
Wash
Elute



Scan for detailed instruction manual



MP BIOMEDICALS

APAC: +65 6775 0008 | custserv.ap@mpbio.com
EUROPE: 00800 7777 9999 | custserv.eur@mpbio.com
AMERICAS: 800 854 0530 | custserv.na@mpbio.com

Learn more at www.mpbio.com

