

\* FOR RESEARCH USE ONLY  
 \* STORE AT 4°C UPON ARRIVAL

## Instruction manual

### Zinc Assay kit LS

#### ( 5-Br-PAPS Chromogenic method )

**Description**

This product is a direct colorimetric assay kit without deproteinization of the sample. At alkaline pH, in a buffered media, zinc reacts with the specific complex, forms a stable colored complex. The color intensity is proportional to the amount of zinc present in the sample.

Zinc is a cofactor of more than 200 kinds of metalloenzymes, and is also a trace element concerned in synthetic of ribonucleic acid or protein. It is widely known that acute zinc deficiency during the growth stage of the mammalian, results in a severe impairment of the skin or hair, and may lead to arrested development. Zinc is essential to reproduction of cell and its relevant supply is necessary for healthy growth.

**Kit contents**

50 tests (Catalog # : ZN01ME )

R-A Buffer ●	12 mL×1
R-R Chelate color (5-Br-PAPS) ●	0.27 mL×1
STD Zinc Standard 200 µg/dL ●	1.65 mL×1

(Catalog # : ZN02ME)=(Catalog # : ZN01ME) ×2

**Note**

- A) Unstableness of incubation temperature may result in unstable results.
- B) Use disposable test tube and glassware washed with 1M HNO<sub>3</sub> or 1M HCl solution and distilled water.
- C) Accuracy in pipetting volume for samples and reagents may affect the quality of assay. Please note that samples, standards and Working Reagent must be poured accurately µL level.
- D) Temperature for chromogen reaction may affect optical density. Please try to extend or shorten chromogen reaction time depending on room temperature.
- E) In the cell lysate or the tissue extract use as specimen, high concentration of proteins or lipid, may affect observed value. Please remove its by ultrafiltration or centrifugation.
- F) Species of zinc-porphyrins cannot be measured in this assay kit.

**Operation**

**1. Sample preparation**

◇Serum or Plasma

Insoluble substances in serum and plasma samples should be removed by filtration or centrifugation. EDTA-plasma cannot be used.

◇Tissue extract, Lysate, Other samples.

Urine (24 hour pooled urine), or other biological fluid:

Add 6M HCl to the sample and adjust pH 2.0-3.0 (e.g. 5-10µL 6M HCl/ 1mL of lysate.). Centrifuge at 6,000 rpm for 15 min. Collect the supernatant and use it for assay.

Tissue:

Add 3% TCA solution, vortex 1 min. and incubate at 4-8°C for 30 min. Centrifuge at 6,000 rpm for 15 min. Collect the supernatant and use it for assay.

\* Sample pH should be between pH2 to pH8.

**2. Assay preparation**

**(1)Bring all reagents to room temperature before use.**

**(2)Prepare enough working Reagent (WR).**

	1 test	Example: 50 tests
R-A Buffer ●	230 (µL)	11.5 (mL)
R-R Chelate color ●	5 (µL)	250 (µL)

\* WR is stored at 2-8 °C and use within one month after prepared.

**3. Assay procedure.**

**Procedure using microplate reader.**

**(1 assay sample 242μL)**

**○ Assay**

- (1) Add 12 μL of Distilled water (Blank) / STD (Standard)/ sample into each well.
- (2) Add 230 μL of Working Reagent ( WR ) to each well and incubate at room temperature for 5 min.
- (3) Read the absorbance at 560 nm (main) and 700nm(sub).  
 --> OD  
 (Sensitivity: 550nm max, 570nm 60%, 580nm 20% or less)

		Assay Sample		
Add	(μL)	Blank OD <sub>Bl</sub>	Standard OD <sub>Std</sub>	Sample OD <sub>S</sub>
1	Distilled water	12	-	-
	STD	-	12	-
	Assay sample	-	-	12
2	WR	230	230	230

↓

Mix and incubate for 5 minutes at room temperature  
 Read the absorbance at 560 nm (main) and 700nm(sub).

**○ Calculations**

$\Delta OD_{Std} = OD_{Std} - OD_{Bl}$

$\Delta OD_S = OD_S - OD_{Bl}$

**Zinc (μg/dL) =  $\Delta OD_S / \Delta OD_{Std} \times 200$**

**Zinc (μM) =  $\Delta OD_S / \Delta OD_{Std} \times 30.6$**

**(Assay example)**

	OD (560nm)	OD (700nm)	OD	ΔOD	Zinc (μg/dL)
Blank	0.062	0.030	0.032	-	-
Standard	0.206	0.031	0.175	0.143	-
Sample	0.117	0.033	0.084	0.052	72.7

**\*Observed 560 nm with 700 nm**

**[OD = OD(560nm) - OD(700nm) ]**

$\Delta OD_{Std} = ( 0.206 - 0.031 ) - ( 0.062 - 0.030 ) = 0.143$

$\Delta OD_S = ( 0.117 - 0.033 ) - ( 0.062 - 0.030 ) = 0.052$

$Zinc_{Sample} (\mu g/dL) = \Delta OD_S / \Delta OD_{Std} \times 200$   
 $= 0.052 / 0.143 \times 200 = 72.7 (\mu g/dL)$

$Zinc_{Sample} (\mu M) = \Delta OD_S / \Delta OD_{Std} \times 30.6$   
 $= 0.052 / 0.143 \times 30.6 = 11.1 (\mu M)$

**\*Observed 560 nm only**

**[OD = OD(560nm)]**

$\Delta OD_{Std} = 0.206 - 0.062 = 0.144$

$\Delta OD_S = 0.117 - 0.062 = 0.055$

$Zinc_{Sample} (\mu g/dL) = \Delta OD_S / \Delta OD_{Std} \times 200$   
 $= 0.055 / 0.144 \times 200 = 76.4 (\mu g/dL)$

$Zinc_{Sample} (\mu M) = \Delta OD_S / \Delta OD_{Std} \times 30.6$   
 $= 0.055 / 0.144 \times 30.6 = 11.7 (\mu M)$

\*In diluted sample of seminal fluid, multiply the result by dilution-factor.

**Performance**

Measuring range 4.0 - 1,000 μg/dL  
 Imprecision Imprecision was evaluated using commercially available quality control serum.

Within run			
	Mean μg/dL	S.D	C.V %
Level 1	68.96	3.05	4.4
Level 2	109.71	2.45	2.2

Interferences No interference by the note of substances were observed.  
 Conjugated bilirubin and unconjugated bilirubin 15 mg/dL  
 Triglyceride 500 mg/dL

**Expiration date and preservation conditions**

Storage conditions: Store at 2-8°C. Don't freeze.  
 Expiration: 1 year from the date of manufacture.  
 After the bottles are opened, the kit should be used in 1 month.

**Reference**

- (1) Makino. T, Saito. M, Horiguchi. D, and Kina. K : A highly sensitive calorimetric determination of serum zinc using water-soluble pyridylazo dye. *Clinical Chimica Acta*, 120, p127-135 (1982).
- (2) Joshua. C, Jia. H, Hirokazu. H, Besim. Ben-N, Bruce. M, Makoto. O: The Therapeutic Effect on Bone Mineral Formation from Biomimetic Zinc Containing Tricalcium Phosphate (ZnTCP) in Zinc-Deficient Osteoporotic Mice. *PLoS One*, 8(8) (2013)

**Manufacturing-and-selling contractor**

Metallogenics Co.,Ltd. Sales Dept.  
 1-14-13, Fujimi, Chuo-ku, Chiba-shi, Chiba 260-0015 Japan  
 TEL : +81-43-227-6767  
 FAX : +81-43-227-6768  
 e-mail : sales@ak-j.com  
 URL: http://metallogenics.co.jp/