

***Pfu* DNA Polymerase, Economy**

02-031 200 U (2.5U/ul), 02-031-5 5 x 200 U (2.5U/ul)

***Pyrococcus furiosus* DNA polymerase (*Pfu* DNA polymerase)** gene was expressed in *E. Coli* in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and 3' →5' exonuclease (proofreading) activity. The MW is 90 kDa, same as that of the natural *Pfu* DNA polymerase.

- *Pfu* DNA polymerase is thermostable and has low error rates.
- It is suitable for PCR and primer extension reactions that require high fidelity synthesis.
- *Pfu* DNA polymerase-generated PCR fragments are blunt-ended.

Applications:

- 1) cloning
- 2) DNA expression
- 3) site-directed mutagenesis

Storage Conditions:

50mM Tris-HCl (pH 8.2), 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.1% Tween20, 0.1% Igepal CA-630, Store at -20°C

Concentration: 2.5 units/ul, where one unit is defined as the amount of enzyme that can incorporate 10 nmols of dNTPs into an acid-insoluble material in 30 minutes at 72°C when activated salmon sperm DNA was used as template/primer.


Quality Assurance: Greater than 95% of protein determined by SDS-PAGE (CBB staining) (Fig.1)
The absence of endonucleases and exonucleases was confirmed.

PCR Test: Good amplification result was obtained in PCR reaction using λDNA as a template (Fig.2).

Reagents Supplied with Enzyme:

10 x Reaction Buffer (*Pfu*): 200mM Tris-HCl (pH 8.8), 100mM KCl, 100mM (NH₄)₂SO₄, 20mM MgSO₄, 1% TritonX-100, 1 mg/ml BSA

<u>General composition of PCR reaction mixture (total 50 ul)</u>	
<i>Pfu</i> DNA polymerase (2.5 units/ul)	0.5 ul
10 x Reaction Buffer (<i>Pfu</i>)	5 ul
2.5mM (each) dNTPs	4 ul
Template	<500ng
Primer 1	0.2~1.0 uM (final conc.)
Primer 2	0.2~1.0 uM (final conc.)
Sterile distilled water	up to 50 ul



Pfu DNA polymerase

PCR condition

98°C 10sec }
55°C 30sec } 30cycles
72°C 10min }
(2min in the case of
2kb DNA)

lane

M : marker
1 : 2 kb
2 : 4 kb
3 : 6 kb
4 : 8 kb

Typical other
BioAcademia supplier




Fig.1 SDS-PAGE of *Pfu* DNA polymerase

Fig.2 Amplification of λ DNA

Related products: # [02-001](#) Taq DNA Polymerase (+dNTPs) # [02-011](#) Taq DNA Polymerase