

## Proteomics Sample Preparation Kit for Cultured Cell (#FMR-001)

**Proteomics Sample Preparation Kit for Cultured Cell** enables to complete sample preparation for proteomics analysis by LC-MS/MS, especially for MRM analysis with mTRAQ® reagent (AB Sciex Pte. Ltd.). This kit contains the necessary reagents for cell lysis, protease digestion, and alkylation.

**For research use only. Not for use in diagnostic procedures.**

**Avoid freeze-thaw cycles.**

### Components (for 10 preps)

- ① Tris HCl (700 µl)
- ② SDS (250 µl)
- ③ Urea (4.24 g)
- ④ Digestion Buffer (0.4 ml)
- ⑤ Protease A (1 vial)
- ⑥ Protease B (2 vial)
- ⑦ Reagent A (10 vials)
- ⑧ Reagent B (10 vials)
- ⑨ Reagent C (1 vial)
- ⑩ Storage Tubes (20 tubes)

### Reagent Preparation

**Cell Lysis Buffer:** Add 5.7 ml of ultrapure water to the Urea bottle (Component ③) just prior to use, and mix thoroughly. Add 700 µl of Tris HCl (Component ①) and then add 250 µl of SDS (Component ②). Mix thoroughly.

**Protease A:** Add 12 µl of ultrapure water to the vial. Store unused solution at -20 °C for 1 month.

**Protease B:** Add 24 µl of ultrapure water to each vial. Store unused solution at -20 °C for 1 month.

**Reagent A, B:** Add 100 µl of ultrapure water just prior to use.

**Reagent C:** Add 1 ml of ultrapure water and dispense 100 µl to each storage tube. The solution is stable for 1 month at -20 °C.

### Storage

①~④, ⑦~⑨: 4 °C, ⑤⑥: -20 °C, ⑩: RT

### Protocols

#### Step A. Cell Lysis

##### [Materials required but not provided]

- Phosphate buffered saline (PBS)
- Trypsin solution
- DMEM medium with 10%FBS
- BSA standard solution (20 mg/ml)
- BCA assay kit
- Coulter counter
- Probe type sonicator (Frequency 20 kHz, Power

Rating 85 W, Probe Size Φ2 mm, or equivalent

##### [General Protocol]

#### Case A. Cell lysis from detached cells

- (1) Wash cultured cells of 10 cm dish with 10 ml of PBS.
- (2) Add 1 to 2 ml of Trypsin solution, then incubate at 37 °C for 5 min.
- (3) Add 8 ml of DMEM medium with 10%FBS to the culture dish.
- (4) Transfer the detached cell suspension to 15 ml tube.
- (5) Centrifuge and discard supernatant.
- (6) Wash the cell pellet with 10ml of PBS.
- (7) Centrifuge and discard supernatant.
- (8) Repeat step 6 and 7 twice.
- (9) Resuspend the cell pellet in PBS (2 to 4 ml).

**NOTE:** You should optimise the cell detachment steps above, depending on the cells you use.

- (10) Count the cell number of 75 µl cell suspension.
- (11) Transfer approx.  $2 \times 10^6$  cells to 1.5 ml tube, centrifuge and discard supernatant.
- (12) Add 200 µl of **Cell Lysis Buffer**.
- (13) Sonicate 30 sec. for twice.
- (14) Add 200 µl of ultrapure water.
- (15) Sonicate 30 sec. for twice.
- (16) If the solution is still viscous, repeat sonication.

#### Case B. Direct lysis from Cell Culture Dish

- (1) Wash cell culture of 10 cm dish with 10 ml of PBS for 3 times.

**NOTE:** Use cell culture dish of  $3 \times 10^6$  to  $1 \times 10^7$  cells.

- (2) Add 500 µl of **Cell Lysis Buffer** and spread whole area of culture dish.
- (3) Add 500 µl of ultrapure water and mix well.
- (4) Use cell scraper and gather the cells to the edge of the dish.
- (5) Sonicate 30 sec. for 3 times or until the solution become low viscous.
- (6) Transfer the lysate to 1.5 ml tube. If the solution is still viscous, repeat sonication.

#### BCA assay

- a. BSA standards preparation: 10-fold dilution of BSA standard solution (20 mg/ml) with the 2-fold diluted **Cell Lysis Buffer**. Prepare BSA

- standards of 0, 0.5, 1, 1.5, 2 mg/ml using 2-fold diluted **Cell Lysis Buffer**.
- Add 10  $\mu$ l of BSA Standards or sample per well.
  - Dispense 200  $\mu$ l of BCA solution per well.
  - Incubate at 37 °C for 30 min.
  - Read at 570 nm. Calculate protein concentrations.

- Centrifuge 13,000 $\times$ g for 2 min.
- Discard supernatant remaining 100  $\mu$ l aliquot.
- Flush centrifuge and discard the remaining liquid.
- Add 28  $\mu$ l of **Digestion Buffer**.
- Incubate at 65 °C for 15 min.
- Leave the solution to ambient temperature.
- Add 28  $\mu$ l of ultrapure water.
- Conduct following BCA assay.

## Step B. Peptide Sample Preparation

### [Materials required but not provided]

- 100% Methanol
- 80% Methanol
- Chloroform
- BSA standard solution (20 mg/ml)
- BCA assay kit
- mTRAQ<sup>®</sup> Reagent (AB SCIEX, #4374771)
- Isopropanol

### [Protocol]

#### Methanol/Chloroform Precipitation

- Transfer 100 to 200  $\mu$ g protein into **Storage Tube**. Fill with 2-fold diluted **Cell Lysis Buffer** to 200  $\mu$ l. You can store the solution at -20 °C

**NOTE:** Be careful not to exceed 200  $\mu$ l. If the protein concentration is low, divide the cell lysate solution into 2 **Storage Tubes** and conduct following step (2) to (17) for each tube, add 14  $\mu$ l of **Digestion Buffer** to the each tube and then combine the solution.

- Add 600  $\mu$ l of 100% Methanol and vortex.
- Add 200  $\mu$ l of Chloroform and vortex.
- Add 400  $\mu$ l of ultrapure water and vortex.

**NOTE:** The solution become clouded after Step (4).

- Incubate on ice for 30 min.
- Centrifuge 13,000 $\times$ g for 5 min.
- Discard the supernatant.
- Add 600  $\mu$ l of cold 100% Methanol and vortex.
- Centrifuge 2,100 $\times$ g for 5 min.

**NOTE:** Use swing bucket rotor centrifuge for this step.

- Centrifuge 13,000 $\times$ g for 2 min.
- Discard supernatant remaining 100  $\mu$ l aliquot.
- Flush centrifuge and discard the remaining liquid.
- Add 1.5 ml of 80% Methanol and vortex.
- Centrifuge 2,100 $\times$ g for 5 min.

**NOTE:** Use swing bucket rotor centrifuge for this step.

- BSA standards preparation: 2-fold dilution of BSA standard solution (20 mg/ml) with 2-fold diluted **Digestion Buffer**. Prepare BSA standards of 0, 1.25, 2.5, 5, 10 mg/ml using 2-fold diluted **Digestion Buffer**.
- Add 2  $\mu$ l of BSA Standards or sample per well.

**NOTE:** Conduct BCA assay in triplicate.

- Dispense 200  $\mu$ l of BCA solution per well.
- Incubate at 37 °C for 30 min.
- Read at 570 nm. Calculate protein concentrations.

#### Protease digestion

- Add 50  $\mu$ l of ultrapure water to the sample.
- Add 1  $\mu$ l of **Protease A**.
- Incubate at 37 °C for 3 hours.
- Add 100  $\mu$ l of ultrapure water and mix by tapping.
- Add 1  $\mu$ l of **Protease B**.
- Incubate at 37 °C for 3 hours.
- Add 1  $\mu$ l of **Protease B**.
- Incubate at 37 °C overnight.
- Add 2.5  $\mu$ l of **Reagent A**.
- Incubate at 37 °C for 30 min, then incubate at RT for a few minutes.
- Add 2.5  $\mu$ l of **Reagent B**.
- Incubate at RT, dark for 30 min.
- Add 1  $\mu$ l of **Reagent C**.
- Incubate at RT for 30 min.
- Lyophilize at -65 °C.

#### Labeling with mTRAQ<sup>®</sup> Reagent

- Add 20  $\mu$ l of ultrapure water.
- Mix by tapping.
- Add 50  $\mu$ l of Isopropanol and vortex.
- Add 20  $\mu$ l of mTRAQ  $\Delta$ 0 (1.0 U). Vortex and spin down.
- Incubate at RT for 90 min.
- Add 150  $\mu$ l of ultrapure water. Incubate at RT for 30 min.
- Incubate at -80 °C for 1 hour.
- Lyophilize at -65 °C for several hours.
- Store at -80 °C until analysis.
- Just before analysis, reconstitute with ultrapure water at 0.2  $\mu$ g/ $\mu$ l concentration.