

Zin-Pro Capture <Mobile Zinc-Responsive Protein Labeling Reagent>

Catalog No. FDV-0013A/ B

Research use only, not for human or animal therapeutic or diagnostic use.

Product Background

Zn²⁺ is tightly bound to proteins and plays essential roles in precise protein folding, and enzymatic function of metalloproteins as cofactor. A simulation suggests about 10% of proteins are estimated as Zn²⁺-binding proteins¹⁾. Mobile Zn²⁺, which is free form or loosely bound to protein surface, is the second class of zinc population and act as a signaling molecule. Recent evidences show an abnormal increase in mobile Zn²⁺ concentration induce pathological events including neuronal degeneration etc. Although fluorescent probes for real time monitoring of mobile Zn²⁺ in live cells have been developed, these probes could not reach comprehensive characterization of Zn²⁺ distribution mechanism and how Zn²⁺ mobiles in cells. To approach zinc homeostasis in cells and to identify proteins involved in zinc distribution on the molecular level, chemical reagents which can be applied for protein experiments are desirable.

Zin-Pro Capture is a chemically labeling reagent of a fluorescein tag to proximal proteins only in the presence of mobile Zn²⁺ ion (Figure 1). Zin-Pro Capture has high selectivity for Zn²⁺ ion, not showing labeling activity with the other abundant metal ions such as Mn²⁺, Fe²⁺/Fe³⁺, Ca²⁺ and Mg²⁺. Its Zn²⁺-responsive labeling activity enable to monitor mobile Zn²⁺ species in the cell on the molecular basis. After labeling Zn²⁺-proximal proteins with fluorescein tag in living cells by Zin-Pro Capture, cells are lysed and the labeled proteins are subsequently enriched by immunoprecipitation with anti-fluorescein antibody. The labeled and purified proteins are characterized by western blotting with antibodies against protein of interest or comprehensively identified by proteomics mass analysis. Observation of intracellular localization of labeled proteins by Zin-Pro Capture is also available in the fixed cells. Miki *et al.* showed that mobile Zn²⁺ in C6 glioma cells was dramatically released to cytosol from metal binding proteins such as metallothionein by oxidative stress and subsequently accumulated in vesicles²⁾. This study revealed with Zin-Pro Capture (AIZin2 in the paper²⁾) that oxidative stress induced Zn-rich vesicles in glioma cells were identified as COPII, COPI and ER-Golgi intermediate compartment (ERGIC) between ER and Golgi apparatus²⁻³⁾.




Figure 1. Principle of zinc-responsive protein labeling reagent, Zin-Pro Capture

In the absence of Zn²⁺ ion, protein labeling activity of Zin-Pro Capture is poor. On the other hand, in the presence of Zn²⁺, Zin-Pro Capture is activated by mobile Zn²⁺ and shows protein labeling activity of fluorescent tag to proximal proteins.




Figure 2. Overview of intracellular labeling by zinc-responsive reagent Zin-Pro Capture

Chemical labeling activity of Zin-Pro Capture is depend on the amount of Zn²⁺ ion in the cell. Highly concentrated Zn²⁺ increases Zin-Pro Capture activity and proteins located on Zn²⁺ rich compartment are highly labeled. On the other hand, in the Zn²⁺ poor region Zin-Pro Capture shows labeling activity with low efficiency and little proteins could be labeled. The labeled proteins can be purified by immunoprecipitation with anti-fluorescein antibody and applied for western blotting and mass spectrometry etc and also observed by microscopy in the fixed cells.

Description

Catalog Number: FDV-0013A/B

Size : 25 µg (FDV-0013A), 25 µg x 3 (FDV-0013B)

Formulation : C₄₆H₄₀N₆O₁₁

Molecular weight : 852.85 g/mol

Chemical structure : See right figure


Solubility : Soluble in DMSO

Warning: Research use only. Not for use in humans.

Fluorescent characteristics: Excitation/ Emission: 495/515 nm

*Compatible with FITC filter

Note: This reagent bears diacetylfluorescein (FDA) which shows little fluorescence emission. FDA promotes cell-permeability of the reagent. FDA is cleaved by intracellular esterases and converted to fluorescein with strong green fluorescence emission.



Application

- SDS-PAGE with fluorescent detection
- Western blotting with antibodies of interest following immunoprecipitation with anti-fluorescein antibody
- Proteomics by MS spectrometry (LC-MS/MS) following immunoprecipitation with anti-fluorescein antibody
- Immunocytochemistry with anti-fluorescein

Reconstitution and Storage

Reconstitution: Stock solution recommended concentration 0.1-1 mM in 100% DMSO.

Storage (powder): Store powder at -20°C.

Storage (solution): After reconstitution in DMSO, aliquot and store at -20°C.

Avoid repeated freeze-thaw cycles and protect from light.

How to use

General procedure for enrichment of labeled proteins

1. Prepare 1 µM of Zin-Pro Capture in the medium and replace medium just before use
2. Culture cells for 20 min with Zin-Pro Capture-containing medium
3. Remove the medium and wash cells with PBS twice
4. Add proper volume of lysis buffer and resuspended lysate well
5. After centrifugation, the supernatants is applied to desalting tools such as dialysis tube, desalting column etc. to remove extra Zin-Pro Capture reagents and their hydrolyzed products.
6. After desalting, add proper amount of anti-fluorescein antibody to lysate and incubate at 4°C
7. Add proper amount of protein A/G agarose beads and incubate at 4°C
8. Wash beads by proper buffers several times
9. Elute the labeled proteins by 50 mM Glycine-HCl buffer etc.

NOTE: Empirically optimize and determine the concentration, incubation time of Zin-Pro Capture and protocol of immunoprecipitation for your experiments.


General procedure for immunocytochemistry

1. Prepare 1 μM of Zin-Pro Capture in the medium and replace medium
2. Culture cells for 20 min with Zin-Pro Capture
3. Remove the medium and wash cells with PBS
4. Cells were fixed with ice-cold MeOH
5. Following step-4, general procedures of immunocytochemistry is available

Application data


Detection of labeled proteins by SDS-PAGE

C6 glioma cells were treated with 200 μM SNOC, a NO generator, for 0, 15 and 180 min. After NO-stimulation, cells were incubated with 1 μM Zin-Pro Capture for 10 min in the absence or presence of 400 μM TPEN, a potent Zn-chelator. Cells were lysed by SDS-sample buffer and proteins were separated in SDS-PAGE and detected by fluorescent imager for labeled proteins or CBB staining for total proteins. Band pattern of labeled proteins was dramatically changed by SNOC treatment and fluorescent signals were suppressed in the presence of TPEN.




Detection of labeled proteins by IP-blot

C6 glioma cells were treated with 200 μM SNOC for 0, 10, and 180 min and incubated with 10 μM Zin-Pro Capture for 10 min. After lysis of cells, labeled proteins were immunopurified by anti-fluorescein antibody (anti-FL). Isolated proteins by anti-FL were separated by SDS-PAGE and detected by western blotting with antibodies of interest (IP-blot). In this case, mitochondrial protein (citrate synthetase), ER protein (calnexin) and Golgi apparatus protein (GRP94) were analyzed. These results indicate mobile Zn was released from metallothionein by NO treatment, transiently accumulated in mitochondria within 10 min and finally transported to ER and Golgi apparatus.



Observation of intracellular localization of labeled protein

C6 glioma cells were treated with 200 μ M SNOC for 10 and 180 min and incubated with 10 μ M Zin-Pro Capture for 10 min. After MeOH-fixation of cells, cells were observed by fluorescent microscopy with the commercial FITC filter set. Localization of labeled proteins which are proximal to high concentrated mobile zinc ion are dramatically changed by treatment time of NO.



Reference

- 1) Andrein *et al.*, *J. Proteome Res.*, **5**, 196–201 (2006) Counting the Zinc-Proteins Encoded in the Human Genome
- 2) Miki *et al.*, *Nat. Methods.*, **13**, 931-937 (2016) A conditional proteomics approach to identify proteins involved in zinc homeostasis
- 3) Motiwala and Martin, *Nat. Methods.*, **13**, 917-918 (2016) Trail-bIAIZin new directions for conditional proteomics

Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter “Product information”) have been prepared based on published research reports on April, 2017. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本製品は学術研究成果を基にフナコシ株式会社が製品化したもので、2017年4月時点における公開研究報告を基に広告文章およびマニュアル(以下、製品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、製品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本製品の仕様および製品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただけますようお願い申し上げます。



E-mail Newsletter
Sign Up

Japanese	English
	

Related products

CellFluor™ GST <Cell-based GST Activity Assay Reagent >


CellFluor™ GST is a novel fluorescent probe for monitoring wide GST members' activity both *in cellulo* or *in vitro*. CellFluor™ GST releases green fluorophore rhodamine 110 upon GST activities. This probe has cell-permeability and can detect intracellular GST activity.

Catalog No. FDV-0030

Size 0.1 μmol

Features

- Easy and quick protocol
- Broad specificity for various GST family members
- Ex/Em: 496 nm/520 nm
(Compatible with commercial FITC filters)



CNBSF <Irreversible GST Inhibitor >

CNBSF is a novel GST inhibitor which irreversibly blocks GST enzymes. CNBSF has membrane-permeability and can be applied into live cell experiments.

Catalog No. FDV-0031

Size 10 mg

Features

- Membrane-permeable and irreversible inhibitor
- Broad specificity for various GST family members
- Covalent inhibition of GSTP1 was experimentally confirmed by MS analysis