

HistoBright <Tissue Structure Fluorescent Dye>

Catalog NO. FDV-0051

Research use only, not for human or animal therapeutic or diagnostic use.

This product has been commercialized under the license from Kochi University and Ehime University.

Product Background

Tissue structures are historically visualized by staining with various colorimetric dyes such as Hematoxylin and Eosin (HE) stains. These conventional colorimetric stains generally require thinly sliced tissue sections and are unsuitable for thick tissue blocks. However, the preparation of thinly sliced tissues loses the 3D structural information. To reconstitute 3D information, a lot of additional thin slices and stainings are required. HistoBright (original compound name; PC in Ref.1) is originally synthesized by Dr. Yosuke Niko, Kochi University, and Dr. Takeshi Imamura, Dr. Masamoto Murakami, and Dr. Ryosuke Kawakami, Ehime University, promotes its biological 3D imaging application. HistoBright is a novel solvatochromic fluorescent dye that changes fluorescent colors in response to solvents or microenvironments and has a high affinity to various biological membranes. The combination of two features of HistoBright contributes to clear and bright imaging of tissue structures.

As HistoBright can be excited by not only one-photon but also two-photon, HistoBright is compatible with two-photon microscopy and useful for tissue imaging. Furthermore, HistoBright, coupled with several tissue-clearing methods such as LUCID and RapiClear, which do not require lipid removal, contributes to analyzing 3D deep tissue imaging of thick tissue slices or blocks, etc., without breaking pathophysiological tissue structures. HistoBright is an innovative fluorescent dye for the visualization of 3D tissue structures.




Figure 1 Overview of HistoBright staining

Description

Catalog Number: FDV-0051

Size: 0.1 mg

Formulation: C₂₇H₂₇NO

Molecular weight: 381.5 g/mol

Solubility: Soluble in DMSO

Fluorescent characteristics:

Ex. 350-500 nm (maximum ~410 nm)

Compatible with 405-488 nm laser (confocal microscopy), 960-1100 nm laser (two-photon microscopy)

Em. 500-800 nm (depending on solvents)

Reconstitution and Storage

Reconstitution: Stock solution recommended concentration 1 mM in 100% DMSO.

Storage (powder): Store powder at -20°C.

Storage (solution): After reconstitution in DMSO, aliquot and store at -20°C.

Avoid repeated freeze-thaw cycles.

How to use

General procedure for tissue staining

*This procedure is an example of HistoBright

1. Tissue samples were excised and then fixed in freshly prepared paraformaldehyde.
2. Tissue blocks or thick slices prepared by any appropriate methods.

NOTE: As HistoBright is a membrane staining dye, any de-lipid procedures including organic solvents or high concentration of detergents such as Triton X100, Tween 20 etc. do not recommended. If you need to use detergents for permeabilization etc., please optimize experimental condition.

3. Stock solution of HistoBright (1 mM) was mixed in tissue clearing solution (LUCID or RapiClear etc.) to prepare the staining and clearing solution (hereafter called Master mix). The initial recommendation of HistoBright concentration is 10 μM.

4. Immerse tissue blocks/sections in the Master mix for over 48 hours at room temperature.


5. After staining tissue blocks/sections, image by confocal laser microscopy or two-photon microscopy.

NOTE: Empirically optimize staining conditions (tissue sample size, pre-treatment of tissues, the concentration of HistoBright, clearing methods, staining time, etc.) for your experiments.

Reference data


Solvatochromic property of HistoBright

Absorbance and fluorescent spectrum of HistoBright in various solvents. While the absorbance spectrum of HistoBright is not sensitive to solvent polarity, HistoBright changes fluorescent color from green to far-red depending on solvent polarity. Under a low polarity environment, HistoBright emits green color. On the other hand, under a high polarity environment, HistoBright emits red to far-red with low quantum yield. Quantum yields are 0.82 (Toluene), 0.77 (EtOAc), 0.83 (CHCl₃), 0.65 (MeCN), 0.20 (EtOH), and 0.07 (MeOH).



Fluorescent spectrum on a liposome membrane model


HistoBright rapidly accumulates in lipid membranes and emits different fluorescent colors in response to the membrane environments. On the liquid order (Lo) model liposome (sphingomyelin (SM)/cholesterol (Chol)), HistoBright emits a maximum of around 575 nm fluorescence with a high quantum yield value (0.72). On the other hand, HistoBright on the liquid disorder (Ld) model liposome (DOPC) exhibits significantly broad and relatively weak fluorescence (quantum yield 0.37), and the maximum wavelength is red-shifted by approximately 30 nm compared to the Lo model.




Application data

Two-photon fluorescence microscopy imaging of human normal skin tissue

After fixation of human skin tissue blocks with 4% paraformaldehyde/PBS, 500 μm thick slides were prepared and treated with the LUCID clearing reagent and 10 μM HistoBright for 76 hours. The tissue slides were observed by two-photon fluorescence microscopy (Ex. 960 nm, Em. cyan channel: 492 nm as SHG, green channel: 500-550 nm, orange channel: 560-593 nm, and red channel: 593-690 nm). The merged image was obtained from four channels.




3D images were reconstructed from Z-stack sequences (step size, 5 μm) of 2D images.




Comparison between confocal microscopy and two-photon microscopy imaging

After fixation of human skin tissue blocks with 4% paraformaldehyde/PBS, 500 μm thick slides were prepared and treated with the RapiClear clearing reagent and 10 μM HistoBright for 72 hours without any permeabilization treatments. Tissues were observed by confocal microscopy (488 nm laser) and two-photon microscopy (960 nm laser).



HE staining after HistoBright fluorescent imaging

After fluorescent imaging, the stained tissue block was washed with PBS. Thin tissue slices were prepared from the tissue block and re-stained with the conventional HE staining procedure.



Notes

All spectrum data were obtained by Dr. Yosuke Niko, Kochi University.

All cell imaging data were obtained by Dr. Ryosuke Kawakami, Ehime University.

Reference

1. Inoue *et al.*, *J. Mater. Chem. B*, 10, 1641-1649 (2022) Synthesis and photophysical properties of a new push-pull pyrene dye with green-to-far-red emission and its application to human cellular and skin tissue imaging.

Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter “Product information”) have been prepared based on published research reports on April, 2023. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本製品は学術研究成果を基にフナコシ株式会社が製品化したもので、2023年4月時点における公開研究報告を基に広告文章およびマニュアル(以下、製品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、製品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本製品の仕様および製品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただけますようお願い申し上げます。

