

Human Mesenchymal Stem Cell (hMSC) Serum-Free Medium for cell proliferation suppression

## Xyltech™ BMT-01

Catalog Number (BBARL/NIPRO): 10201/87-283                    100 mL

### 1. Product features

This product is a synthetic culture medium suitable for suppressive growth control of human mesenchymal stem cell (hMSC). Xyltech™ BMT-01 can be used to control the growth rate of hMSC in combination with Xyltech™ Growth BMT, which is a serum-free medium for hMSC proliferation. This product is serum-free culture medium.

### 2. Precautions for use

Xyltech™ BMT-01 does not contain substances that neutralize trypsin activity. When subculturing cells with trypsin, it is strongly recommended that the trypsin activity be sufficiently neutralized with a trypsin inhibitor. Dilution washing alone does not completely remove trypsin activity and the remaining protease activity will reduce subsequent cell growth. This product is a research reagent. It cannot be used for human or animal treatment or diagnostic purposes.

### 3. Storage

Store Xyltech™ BMT-01 in a cool, dark place (2-8°C). Do not freeze the medium to avoid deterioration of some active ingredients.

### 4. Example of suppressive growth control protocol for normal human adipose-derived stem cells (ADSCs) culture using Xyltech™ BMT-01

#### 4-1. Cells and reagents

- |  |                                      |
|--|--------------------------------------|
| • Normal human ADSCs (100 mm-dish, Fibronectin coated) |                                      |
| • Xyltech™ BMT-01 (growth suppression medium)          | *This product                        |
| • Xyltech™ Growth BMT (growth medium)                  | (CAT. # (BBARL/NIPRO): 10211/87-286) |
| • r-TE (r-Trypsin/EDTA Solution)                       | (NIPRO CAT. #: 87-974)               |
| • s-TI (Synthetic Trypsin Inhibitor Solution)          | (NIPRO CAT. #: 87-975)               |
| • D-PBS (-)  |                                      |

#### 4-2. Growth control of normal human ADSCs

1. Warm the culture medium, D-PBS (-), r-TE, and s-TI in a 37°C water bath.
2. Remove the culture supernatant of selected normal human ADSCs that have reached around 80% confluence (sub confluence).
3. Rinse the cell layer with 5 mL of D-PBS (-).
4. Add 0.5 mL of r-TE and incubate at 37°C for approximately 2 minutes.
5. Add 0.5 mL of s-TI, mix well, gently pipette up and down several times, collect cells from the dish,

- and centrifuge at 1,000 rpm, for 5 minutes.
6. Aspirate the supernatant and add the appropriate amount of Xyltech™ Growth BMT medium to resuspend the cells and seed into a fibronectin coated new tissue culture dish.
  7. The next day, use a phase-contrast microscope to confirm the cells are engrafted and replace the growth medium with the Xyltech™ BMT-01 culture medium.
  8. The cells can be cultured up to 3 days with Xyltech™ BMT-01. The cells start to regrow quickly and become confluent 1-2 days after changing the medium back to the Xyltech™ Growth BMT medium. Begin subculture and/or experiments with the cells.

#### **4-3. Phase contrast microscope images of normal human ADSCs cultured with Xyltech™ Growth BMT or Xyltech™ BMT-01**




Fig. 1 Normal human ADSCs cultured with Xyltech™ Growth BMT for 3 days.




Fig. 2 Normal human ADSCs cultured with Xyltech™ BMT-01 (growth suppressive medium) for 3 days.

Bars=200 μm

\*The protocol is based on experimental results. It may be necessary to adjust seeding density, and passage timing according to the cells. This protocol is intended for research purposes only.

#### **5. For Inquiries about products**

Bourbon Biomedical Advanced Research Laboratories, Inc. (BBARL, Inc.)  
1-3-1, Ekimae, Kashiwazaki-city, Niigata,  
TEL: +81-257-23-2769 E-mail: support@bourbon-barl.co.jp