

Platinum-A Retroviral Packaging Cell Line, Amphotropic

CATALOG NUMBER: RV-102

STORAGE: Liquid nitrogen

Note: For best results begin culture of cells immediately upon receipt. If this is not possible, store at -80°C until first culture. Store subsequent cultured cells long term in liquid nitrogen.

QUANTITY & CONCENTRATION: 1.0 mL, $>3 \times 10^6$ cells/mL in DMEM, 20% FCS and 10% DMSO

Background

Retroviruses are efficient tools for delivering heritable genes into the genome of dividing cells. However, conventional NIH-3T3 based retroviral packaging cell lines have limited stability and produce low viral yields, mainly due to poor expression level of the retroviral structure proteins (gag, pol, env) in the packaging cells.

The Platinum-A (Plat-A) Cell Line, a potent retrovirus packaging cell line based on the 293T cell line, was generated using novel packaging constructs with an EF1 α promoter to ensure longer stability and high-yield retroviral structure protein expression (gag, pol, amphotropic env). Plat-A cells can be kept in good condition in for at least 4 months in the presence of drug selection, and can produce retroviruses with an average titer of 1×10^6 infectious units/mL by transient transfection. In addition, replication competent retroviruses (RCR) are virtually nonexistent because only coding sequences of viral structural genes are used, avoiding any unnecessary retroviral sequences.

The Plat-A cell line is designed for rapid, transient production of high-titer, amphotropic retrovirus.




Figure 1. TF-1 cells were infected with GFP retrovirus supernatant produced in Plat-A cells after transfection with pMX-GFP.

Medium

1. Culture Medium: DMEM, 10% fetal calf serum (FCS), 1 µg/mL puromycin, 10 µg/mL blasticidin, penicillin and streptomycin
2. Freeze Medium: 70% DMEM, 20% FBS, 10% DMSO

Methods

I. Establishing Plat-A Cultures from Frozen Cells

1. After quickly thawing the cells in a 37°C water bath, immediately transfer the thawed cell suspension into a 15 mL tube containing 10 mL of culture medium.
2. Centrifuge the tube for 5 min at 1300 to 1500 rpm.
3. Discard the supernatant and break the cell pellet by finger tapping.
4. Add a few drops of culture medium with gentle shaking and finger tap the tube a few times.
5. Add 2 mL of culture medium to the tube and gently pipet the cell suspension up and down twice.
6. Transfer the cell suspension to a 10 cm culture dish (Falcon® #3003 works well) containing 8 mL of culture medium.
7. Swirl the culture plate well to mix the cells, then incubate the cells for three days before expansion.

Important Notes:

- ***Don't change the culture medium during the first three days. It is normal to see some cells floating after the first 24 hours.***
- ***Don't culture cells to complete confluency. Split cells 4X to 6X every two to three days when the culture reaches 70-90% confluency.***

II. Splitting the Cells

Note: Avoid forming bubbles as much as possible during this procedure.

1. Wash cells once with PBS.
2. Add 4 mL of 0.05% Trypsin/0.5 mM EDTA solution to a 10 cm dish and incubate at 37°C for 3-5 min.
3. Remove the cells from the dish surface by tapping the rim of the culture dish.
4. Transfer 10 mL of the culture medium to a 50 mL tube.
5. Using the same pipette with some residual culture medium, wash the dish surface gently three times in 4 mL of the Trypsin/EDTA solution.
6. Gently pipette the cell suspension up and down 7 times and transfer the cell suspension into the 50 mL tube containing 10 mL medium from step 4.
7. Centrifuge the cells for 5 min at 1300-1500 rpm.
8. Discard the supernatant and break the cell pellet by finger tapping.
9. Add a few drops of culture medium with gentle shaking and finger tap the tube a few times.
10. Add 5 mL of culture medium and gently pipet the cell suspension up and down twice.
11. Add 15 mL of culture medium, then count and seed the cells. Typically 10^7 cells can be harvested from one 10 cm culture dish.

Transfection

1. Seed 2×10^6 cells in a 60 mm culture dish without antibiotics including puromycin and blasticidin one day before transfection.
2. After 16 to 24 hours, start transfection when the culture becomes 70-80% confluent.
Note: We suggest transfecting cells with FuGENE® Transfection Reagent (Roche Applied Science) or Lipofectamine™ Plus (Invitrogen). For example, 3 µg retroviral expression plasmid is mixed with 9 µL FuGENE® Transfection Reagent according to the manufacturer's recommendation. The mixed DNA- FuGENE® complex is added by dropwise into the culture media.
3. Harvest retroviral supernatant 48 hours after transfection.

References

1. Morita, S., Kojim, T., and Kitamura, T. (2000) *Gene Therapy* 7: 1063-1066.

Recent Product Citations

1. Lee, J. H. et al. (2015). AKT phosphorylates H3-threonine 45 to facilitate termination of gene transcription in response to DNA damage. *Nucleic Acids Res.* doi: 10.1093/nar/gkv176.
2. Han, F. et al. (2015). Human induced pluripotent stem cell–derived neurons improve motor asymmetry in a 6-hydroxydopamine–induced rat model of Parkinson's disease. *Cytotherapy.* doi: 10.1016/j.jcyt.2015.02.001.
3. Balboni, A. L. et al. (2015). p53 and DeltaNp63alpha co-regulate the transcriptional and cellular response to TGFbeta and BMP Signals. *Mol Cancer Res.* doi: 10.1158/1541-7786.MCR-14-0152-T.
4. Ji, Y. et al. (2015). miR-155 augments CD8+ t-cell antitumor activity in lymphoreplete hosts by enhancing responsiveness to homeostatic γ C Cytokines. *Proc Natl Acad Sci U S A.* **112**:476-481.
5. DeCastro, A. J. et al. (2015). Δ NP63 α transcriptionally activates chemokine receptor 4 (CXCR4) expression to regulate breast cancer stem cell activity and chemotaxis. *Mol Cancer Ther.* **14**:225-235.
6. Hu, R. et al. (2015). NF- κ B signaling is required for XBP1 (unspliced and spliced)-mediated effects on antiestrogen responsiveness and cell fate decisions in breast cancer. *Mol Cell Biol.* **35**:379-390.
7. Robin, J. D. et al. (2015). Isolation and immortalization of patient-derived cell lines from muscle biopsy for disease modeling. *J Vis Exp.* **95**:e52307-e52307.
8. Hamabe, A. et al. (2014). Role of pyruvate kinase M2 in transcriptional regulation leading to epithelial–mesenchymal transition. *Proc Natl Acad Sci U S A.* **111**:15526-15531.
9. Krivega, I. et al. (2014). Role of LDB1 in the transition from chromatin looping to transcription activation. *Genes Dev.* **28**:1279-1290.
10. Nam, Y.J. et al. (2013). Reprogramming of human fibroblasts toward a cardiac fate. *PNAS.* **110**:5588-5593.
11. Hrdlickova, R. et al. (2012). Alternatively spliced telomerase reverse transcriptase variants lacking telomerase activity stimulate cell proliferation. *Mol. Cell. Biol.* **32**:4283-4296.
12. Wright, L.N. et al. (2012). Modeling the transcriptional consequences of epidermal growth factor receptor ablation in Ras-initiated squamous cancer. *Clin. Cancer Res.* **18**:170-183.
13. Qiao, Y. et al. (2011). FOXQ1 regulates epithelial-mesenchymal transition in human cancers. *Cancer Res.* **71**:3076-3086.
14. Hudecek, M. et al. (2010). The B-cell tumor-associated antigen ROR1 can be targeted with T-cells modified to express a ROR1-specific chimeric antigen receptor. *Blood.* **116**:4532-4541.

15. Oda, Y. et al. (2010). Induction of pluripotent stem cells from human third molar mesenchymal stromal cells. *J. Biol. Chem.* **285**:29270-29278.
16. Wu, Y. et al. (2009). Improved coinfection with amphotropic pseudotyped retroviral vectors. *J. Biochem. Biotechnol.* 10.1155/2209-901079.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS 's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

License Information

This licensed product is intended for ACADEMIC, GOVERNMENT AND NON-PROFIT RESEARCH USE ONLY; not for use in diagnostic or therapeutic procedures. The product may not be transferred, sold, or otherwise provided to another laboratory except by an authorized distributor of Cell Biolabs, Inc.

Use of this product by Biotechnology and Pharmaceutical companies requires a license for all fields of use including research. Please contact:

Director of Business Development
Cell Biolabs, Inc.
busdev@cellbiolabs.com

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2008-2015: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.