

### **Reagents Needed**

Component	Vendor	Catalog #
Calbryte™ 520 AM	AAT Bioquest	20651
PowerLoad™ Concentrate, 100X	Thermo Fisher	P10020
DMEM/F12 (Imaging Medium)	Gibco™	11330057

This protocol is an adapted version of AAT Bio's Calbryte<sup>™</sup> 520 AM product manual which can be found <u>here</u>.

## Resuspending Calbryte 520 AM:

- **1**. Briefly centrifuge 50  $\mu$ g tube
- Resuspend in 18 uL DMSO to make a ~2.5 mM (500X) stock solution
- 3. Store leftover solution at -20°C

#### **Preparing 2X Loading Solution:**

- 1. Calculate the appropriate amount of 2X Loading Solution needed
- 2. Add 1:250 Calbryte<sup>™</sup> 520 AM
- 3. Add 1:50 PowerLoad<sup>™</sup> Concentrate

Notes: 1) The resulting concentration of Pluronic F-127 in PowerLoad<sup>™</sup> Concentrate is unknown. 2) We <u>do not</u> use Probenecid due to potentially toxic effects on the cell and the Calbryte protocol states that the calcium indicator was specifically designed to eliminate the need for probenecid.

#### Vortex

Add appropriate amount of DMEM/F12

#### Loading the Dye:

- 1. Remove half the existing medium from the well. For 96-well plates, leave 50 uL of medium remaining.
- 2. Add the same amount removed of 2X Loading Solution. For 96well plates, add 50 uL of Loading Solution.

Note: The final concentration of Calbryte<sup>™</sup> 520 AM should be 5 uM and PowerLoad<sup>™</sup> Concentrate 1X.

3. Incubate for 40-50 minutes at 37 °C.

Note: Incubation for 30 and 70 minutes also works, but signal and background may de/increase.

- Exchange the Loading Solution with DMEM/F12 (imaging medium) twice. For 96-well plates, add 90 uL of DMEM/F12.
- 5. Incubate the cells for 15 minutes at room temperature.

Note: This "incubation" time typically occurs while setting up everything at the microscope, finding a sweet spot of cells, adjusting focus, etc.

# Calcium Imaging Quick Guide RealDRG<sup>™</sup> Nociceptors

#### Imaging Setup

Fluorescent videos of RealDRG<sup>™</sup> cultures were acquired using a Leica DMI6000B microscope, a DFC365FX camera, and a 10x air objective (NA 0.25, Leica Microsystems) for a duration of 185 seconds at 4 frames per second.

### Drug Dosing:

- 1. Drugs were prepared 10X working concentration in 1% DMSO in DMEM/F12.
- 2. Image acquisition was started on a random field of view with >200 cells.
- 3. Within 10 seconds of starting image acquisition, drugs were added to the cultures. For 96-well plates, 10 uL of the drug was added.

Comp	ound	List:
------	------	-------

Drug	Working Conc.	Vendor	Catalog #
Veratridine	1 µM	Tocris	2918
Tetrodotoxin	1 µM	Biotium	00061
KCI	30 mM	Millipore Sigma	P5405
Capsaicin*	100/500 nM	Sigma-Aldrich	M2028
Capsazepine*	1 µM	Tocris	0464
Menthol*	100 µM	Sigma-Aldrich	M2772

Note: RealDRG<sup>™</sup> Nociceptors only respond to the compounds marked with an asterisk after a minimum of 4 weeks of maturation in Senso-MM.

#### Image Analysis:

Analysis was automated using FIJI

- 1. Individual somata were detected by local maxima-based image segmentation and binary threshold
- 2. Mean intensity was measured over time inside each soma
- 3. The intensity profiles were graphed, in terms of  $\Delta$ F/F0, where  $\Delta$ F is the intensity differential compared to minimum fluorescence (F0) of the given sample