Challenge

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Operating the CL-3000 as part of an approved UVGI Workflow

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We do not advocate specific treatments or approaches. We are simply sharing the most recent evidence from the medical community to help HCWs during the SARS-CoV-2 pandemic. Your UVGI workflow should be set and approved by your hospital/institution. Users can refer to the UNMC² and/or Stanford University³ workflows for process recommendations.

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References

1. CDC. Coronavirus Disease 2019 (COVID-19) - Decontamination and Reuse of filtering Facepiece Respirators. Centers for Disease Control and Prevention <https://www.cdc.gov/coronavirus/2019-ncov/hcp/ppe-strategy/decontamination-reuse-respirators.html> (2020).
2. CLowe, J. et al. N95 Filtering Facemask Respirator Ultraviolet Germicidal Irradiation (UVGI) Process for Decontamination and Reuse. <https://www.nebraskamed.com/sites/default/files/documents/covid-19/n-95-decon-process.pdf> (2020).
3. Price, A. & Chu, L. COVID-19 Evidence Service - Addressing COVI-19 Face Mask Shortages [v1.2]. <https://aim.stanford.edu/covid-19-evidence-service/> (2020).
4. Darnell, M. E. R., Subbarao, K., Feinstone, S. M. & Taylor, D. R. Inactivation of the coronavirus that induces severe acute respiratory syndrome, SARS-CoV. *Journal of Virological Methods* 121, 85–91 (2004).
5. McDevitt, J. J., Milton, D. K., Rudnick, S. N. & First, M. W. Inactivation of Poxviruses by Upper-Room UVC Light in a Simulated Hospital Room Environment. *PLoS ONE* 3, e3186 (2008).
6. Rutala, W. A., Gergen, M. F. & Weber, D. J. Room Decontamination with UV Radiation. *Infect. Control Hosp. Epidemiol.* 31, 1025–1029 (2010).
7. Fisher, E. M. & Shaffer, R. E. A method to determine the available UV-C dose for the decontamination of filtering facepiece respirators. *Journal of Applied Microbiology* 110, 287–295 (2011).
8. Park, G. w., Linden, K. g. & Sobsey, M. d. Inactivation of murine norovirus, feline calicivirus and echovirus 12 as surrogates for human norovirus (NoV) and coliphage (F+) MS2 by ultraviolet light (254nm) and the effect of cell association on UV inactivation. *Letters in Applied Microbiology* 52, 162–167 (2011).
9. Zou, S. et al. Inactivation of the novel avian influenza A (H7N9) virus under physical conditions or chemical agents treatment. *Virol J* 10, 289 (2013).
10. Kiryu, I., Sakai, T., Kurita, J. & Iida, T. Virucidal Effect of Disinfectants on Spring Viremia of Carp Virus. *Fish Pathol.* 42, 111–113 (2007).
11. Singh, R., Sharma, A., Hong, S. & Jang, J. Electrical immunosensor based on dielectrophoretically-deposited carbon nanotubes for detection of influenza virus H1N1. *Analyst* 139, 5415–5421 (2014).
12. Beck, S. E., Wright, H. B., Hargy, T. M., Larason, T. C. & Linden, K. G. Action spectra for validation of pathogen disinfection in medium-pressure ultraviolet (UV) systems. *Water Research* 70, 27–37 (2015).
13. Lindsley, W. G. et al. Effects of Ultraviolet Germicidal Irradiation (UVGI) on N95 Respirator Filtration Performance and Structural Integrity. *Journal of Occupational and Environmental Hygiene* 12, 509–517 (2015).
14. Bae, K. S., Shin, G.-A., Bae, K. S. & Shin, G.-A. Inactivation of various bacteriophages by different ultraviolet technologies: Development of a reliable virus indicator system for water reuse. *Environmental Engineering Research* 21, 350–354 (2016).
15. Bedell, K., Buchaklian, A. & Perlman, S. Efficacy of an automated multi-emitter whole room UV-C disinfection system against Coronaviruses MHV and MERS-CoV. *Infect Control Hosp Epidemiol* 37, 598–599 (2016).
16. Chu, H.-A. & Chiu, Y.-F. The Roles of Cytochrome b559 in Assembly and Photoprotection of Photosystem II Revealed by Site-Directed Mutagenesis Studies. *Front. Plant Sci.* 6, (2016).
17. Song, K., Mohseni, M. & Taghipour, F. Application of ultraviolet light-emitting diodes (UV-LEDs) for water disinfection: A review. *Water Research* 94, 341–349 (2016).
18. Ziegler, C. M. et al. The Lymphocytic Choriomeningitis Virus Matrix Protein PPLY Late Domain Drives the Production of Defective Interfering Particles. *PLoS Pathog* 12, (2016).
19. Berger, A. K. et al. Viral RNA at Two Stages of Reovirus Infection Is Required for the Induction of Necroptosis. *J. Virol.* 91, e02404-16, /jvi/91/6/e02404-16.atom (2017).

20. Fryk, J. J. et al. Reduction of Zika virus infectivity in platelet concentrates after treatment with ultraviolet C light and in plasma after treatment with methylene blue and visible light. *Transfusion* 57, 2677–2682 (2017).
21. Fung, T. S. & Liu, D. X. Activation of the c-Jun NH 2 -terminal kinase pathway by coronavirus infectious bronchitis virus promotes apoptosis independently of c-Jun. *Cell Death & Disease* 8, 1–13 (2017).
22. Kim, D.-K., Kim, S.-J. & Kang, D.-H. Inactivation modeling of human enteric virus surrogates, MS2, Q β , and Φ X174, in water using UVC-LEDs, a novel disinfecting system. *Food Research International* 91, 115–123 (2017)..
23. Bodmer, B. S., Fiedler, A. H., Hanauer, J. R. H., Prüfer, S. & Mühlebach, M. D. Live-attenuated bivalent measles virus-derived vaccines targeting Middle East respiratory syndrome coronavirus induce robust and multifunctional T cell responses against both viruses in an appropriate mouse model. *Virology* 521, 99–107 (2018).
24. Frascaroli, G. et al. Human Macrophages Escape Inhibition of Major Histocompatibility Complex-Dependent Antigen Presentation by Cytomegalovirus and Drive Proliferation and Activation of Memory CD4+ and CD8+ T Cells. *Front. Immunol.* 9, (2018).
25. Lee, M. K. et al. Transmembrane Protein pUL50 of Human Cytomegalovirus Inhibits ISGylation by Downregulating UBE1L. *J Virol* 92, e00462-18, /jvi/92/15/e00462-18.atom (2018).
26. Mathew, A. M., Mun, A. B. & Balakrishnan, A. Ultraviolet Inactivation of Chikungunya Virus. *Intervirology* 61, 36–41 (2018).
27. Mills, D., Harnish, D. A., Lawrence, C., Sandoval-Powers, M. & Heimbuch, B. K. Ultraviolet germicidal irradiation of influenza-contaminated N95 filtering facepiece respirators. *American Journal of Infection Control* 46, e49–e55 (2018).
28. Nishisaka-Nonaka, R. et al. Irradiation by ultraviolet light-emitting diodes inactivates influenza A viruses by inhibiting replication and transcription of viral RNA in host cells. *Journal of Photochemistry and Photobiology B: Biology* 189, 193–200 (2018).
29. Vaidya, V. et al. Ultraviolet-C irradiation for inactivation of viruses in foetal bovine serum. *Vaccine* 36, 4215–4221 (2018).
30. Yang, L. et al. Porcine Epidemic Diarrhea Virus-Induced Epidermal Growth Factor Receptor Activation Impairs the Antiviral Activity of Type I Interferon. *J Virol* 92, e02095-17, /jvi/92/8/e02095-17.atom (2018)
31. Baidaliuk, A. et al. Cell-Fusing Agent Virus Reduces Arbovirus Dissemination in *Aedes aegypti* Mosquitoes In Vivo. *J Virol* 93, e00705-19, /jvi/93/18/JVI.00705-19.atom (2019).
32. Campbell, T. M. et al. Functional paralysis of human natural killer cells by alphaherpesviruses. *PLoS Pathog* 15, (2019).
33. DeFord, D. M. et al. Evaluation of the role of respiratory syncytial virus surface glycoproteins F and G on viral stability and replication: implications for future vaccine design. *Journal of General Virology*, 100, 1112–1122 (2019).
34. Heimbuch, B. & Harnish, D. Research to mitigate a shortage of respiratory protection devices during public health emergencies. (2019).
35. Lee, A. C. Y. et al. H7N9 influenza A virus activation of necroptosis in human monocytes links innate and adaptive immune responses. *Cell Death & Disease* 10, 1–16 (2019).
36. Romero, N., Van Waesberghe, C. & Favoreel, H. W. Pseudorabies virus infection of epithelial cells leads to persistent but aberrant activation of the NF- κ B pathway, inhibiting hallmark NF- κ B-induced pro-inflammatory gene expression. *J Virol* JVI.00196-20, jvi;JVI.00196-20v1 (2020) doi:10.1128/JVI.00196-20.
37. Zhao, X. et al. Activation of C-Type Lectin Receptor and (RIG)-I-Like Receptors Contributes to Proinflammatory Response in Middle East Respiratory Syndrome Coronavirus-Infected Macrophages. *J Infect Dis* 221, 647–659 (2020).
38. Zaffina, S. et al. Accidental exposure to UV radiation produced by germicidal lamp: case report and risk assessment. *Photochem. Photobiol.* 88, 1001–1004 (2012).
39. International Commission on Non-Ionizing Radiation Protection. ICNIRP statement-Protection of workers against ultraviolet radiation. *Health Phys* 99, 66–87 (2010).

This document is true and correct at the time of publication; the information within is subject to change. Other documents may supersede this document, including technical modifications and corrections.

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© Analytik Jena US LLC | Pictures ©: Analytik Jena / freepik
Printout and further use permitted with reference to the source

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