

Datasheet: SFL1267F

Description:	MOUSE ANTI HUMAN CD4:FITC
Specificity:	CD4
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	RPA-T4
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

PDF Available [Download PDF...](#)

Target Species Human

Product Form See PDF datasheet

Note To view the datasheet in alternative languages, please follow the link below:
[Alternative languages](#)

Please contact the AbD Serotec Technical Support team (tech.uk@abdserotec.com) regarding any questions or comments upon the performance of this product.

Buffer Solution Phosphate buffered saline

Preservative 0.09% Sodium Azide (NaN₃)
Stabilisers 1% Bovine Serum Albumin

External Database Links

UniProt:
[P01730](#) CD4_HUMAN [Related reagents](#)

Entrez Gene:
[920](#) CD4 [Related reagents](#)

Storage Keep in a dark place at 2-8oC or if aliquotted, at -20oC for long term storage. Avoid multiple freeze/thaw cycles.

Shelf Life Please refer to vial label for expiry date

Health And Safety Information Material Safety Datasheet Documentation #10041 available at:
<http://www.abdserotec.com/uploads/MSDS/10041.pdf>

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

Other formats and applications

This reagent is also available in a [RPE](#) format

For research purposes only, unless otherwise specified in writing by AbD Serotec.

'M189831:110310'



AbD Serotec offices contact details click [here](#)

Address:

MorphoSys UK Ltd, Endeavour House, Langford
Business Park, Langford Lane, Kidlington,
Oxford, OX5 1GE, UK

Tel: +44 (0)1865 852700

Email: tech.uk@abdserotec.com

Web: www.abdserotec.com

Printed on 09 Nov 2011