



DATA SHEET

Copyright GeneTex, Inc. All rights reserved.

Endo/Blocker

Catalog Number	GTX73347
Product Description	Endo / Blocker (Endogenous Enzyme Blocker)
Background	<p>This Endogenous Enzyme Blocker is useful to suppress endogenous phosphatase and peroxidase-like enzyme activities found in many tissue preparations destined for immunocytochemical procedures. Routine treatment of tissue specimens with Endo/Blocker aids in the elimination of undesirable background staining when using these immunoenzymatic procedures. The treatment of tissue sections with Endo/Blocker does not result in any noticeable decrease of the antigenic characteristics of the specimen. This endogenous enzyme blocker is useful for suppressing endogenous phosphatase and peroxidase-like enzyme activities found in many tissue preparations destined for immunocytochemical procedures. Routine treatment of tissue specimens with Endo/Blocker aids in the elimination of undesirable background staining when using these immunoenzymatic procedures. The treatment of tissue sections with Endo/Blocker does not result in any noticeable decrease of the antigenic characteristics of the specimen.</p>
Tested Applications	Immunohistochemistry. The usefulness of this product in other applications has not been determined.
Application Note	For IHC: Please refer to protocol information. Not tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
Storage Instruction	Keep as concentrated solution, aliquot and store at 4 [°] C. Do not freeze.
