

BNP-32 (Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-011-10	1 kit	\$420.00	<input type="button" value="ORDER"/>

Product description **More information**

Brain Natriuretic Peptides (BNP)
 Normal Plasma Peptide/Protein Levels

Protocol [Assay Protocol \(PDF\)](#)

Sample Preparation [Blood collection](#)
 [Plasma extraction](#)
 [Tissue extraction](#)
 [CSF extraction](#)

Catalog # EK-011-10

Standard Size 96 wells

Sensitivity Minimum Detectable Concentration = 0.45 ng/ml

Range 0-100 ng/ml


Linear Range 0.45-6.3 ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	BNP-32 (Porcine)	100
	BNP-26 (Porcine)	100
	BNP-32 (Human)	0
	BNP-32 (Rat)	0
	BNP-45 (Rat)	0
	α-ANP (1-28) (Human, Canine)	0

Sequence Ser - Pro - Lys - Thr - Met - Arg - Asp - Ser - Gly - Cys - Phe - Gly - Arg - Arg - Leu - Asp - Arg - Ile - Gly - Ser - Leu - Ser - Gly - Leu - Gly - Cys - Asn - Val - Leu - Arg - Arg - Tyr

Disulfide Bridge: Cys10 - Cys26

Image



Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.