

Dynorphin A (Human, Rat, Mouse, Porcine) - EIA Kit

Catalog#	Standard size	Price	Order
EK-021-03	1 kit	\$420.00	<input type="button" value="ORDER"/>
Product description		More information	References

Beta-Neoeendorphin-Dynorphin Precursor Schematic (Human)


Protocol	Assay Protocol (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	EK-021-03
Standard Size	96 wells
Sensitivity	Minimum Detectable Concentration = 0.11ng/ml
Range	0-100 ng/ml
Linear Range	0.11-0.74ng/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Dynorphin A (Porcine)	100
	Dynorphin A (1-13) (Porcine)	<0.01
	Dynorphin A (1-8) (Porcine)	0
	Dynorphin B (Porcine)	0
	β-Endorphin (Human)	0
	α-Neo-Endorphin (Porcine)	0
	Leu-Enkephalin	0

Sequence Tyr - Gly - Gly - Phe - Leu - Arg - Arg - Ile - Arg - Pro - Lys - Leu - Lys - Trp - Asp - Asn - Gln

Image

Dynorphin A (Pro-Dynorphin (209-225)) (Porcine)



4-P Fit: $y = (A - D)/(1 + (x/C)^B) + D$

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.