

CNP-22 (Human, Rat, Mouse, Porcine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-012-03	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

C-type Natriuretic Peptides

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-012-03

Standard Size 96 well plate


Sensitivity Minimum Detectable Concentration = 6.6 pg/ml

Range 0-10000 pg/ml

Linear Range 6.6-176 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CNP-22 (Human, Rat, Porcine)	100
	CNP-53 (Rat, Porcine)	100
	Urodilatin	8.47
	CNP-53 (1-29) (Rat, Porcine)	0
	Alpha-ANP (1-28) (Human, Canine, Porcine, Bovine)	0
	ANF (Rat, Rabbit, Mouse)	0
	BNP-32 (Porcine)	0
	BNP-45 (Rat)	0
	Pro-BNP (1-46) (Human)	0
	Pro-BNP (22-46) (Human)	0
	Pro-BNP (47-76) (Human)	0
	N-Terminal Prepro-BNP (1-76)(Human)	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Substance P	0
	DNP	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$