

VIP (Human, Rat, Mouse, Porcine, Ovine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-064-16	1 kit	\$460.00	ORDER

Product description

More information

References

[VIP \(Vasoactive Intestinal Peptide /PHI and Related Peptides\)](#)

[Prepro-VIP/PHM-27 Schematic \(Human\)](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-064-16

Standard Size 96 well plate


Sensitivity Minimum Detectable Concentration = 13.3 pg/mL

Range 0-10000 pg/ml

Linear Range 13.3-301 pg/ml

Cross Reactivity	Peptide	Cross-reactivity
	VIP (Human, Porcine, Rat)	100
	VIP (Guinea Pig)	100
	VIP (10-28) (Human, Porcine, Rat)	34.6
	VIP (Chicken)	28
	PACAP-27-NH ₂ (Human, Rat, Ovine)	<0.02
	VIP (1-12) (Human, Porcine, Rat)	0
	PHM-27 (Human)	0
	Substance P	0
	Endothelin-1 (Human, Rat, Porcine, Mouse, Canine, Bovine)	0
	Secretin (Porcine)	0
	Glucagon (Human)	0
	Galanin (Rat)	0
	Somatostatin	0
	PACAP-38 (Human, Rat, Ovine)	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$