

## PRODUCT INFORMATION

---

<b>Product Name :</b>	Sonicated Salmon Sperm DNA		
<b>Code No. :</b>	F012, F013		
<b>Size :</b>	1 ml / tube	F012	1 ml F013 5 ml ( 1ml × 5 )
<b>Concentration:</b>	10 mg/ml		
<b>Description:</b>	Salmon Sperm DNA was shared by sonication. DNA is dissolved in deionized distilled water, with a fragment size range of 300-700 base pairs.		
<b>Quality Control:</b>	Intrinsic DNase activity was not detected in the product, by the use of pUC18 plasmid DNA as a control.		
<b>Preparation:</b>	Boiling for 3 minutes and cooling in an ice bath before use.		
<b>Storage Condition:</b>	-20°C		

**For Research Use Only. Not for use in diagnostic procedures.**

Ver. 2.1