

## RISH™ Retrieval, 10X

Pretreatment Reagent

Control Number: 902-0209-090717

**Catalog Number:** BRR 0209 M

**Description:** 500 ml; concentrate

**Intended Use:**

For Research Use Only. Not for use in diagnostic procedures.

**Summary & Explanation:**

RISH™ Retrieval is a heat retrieval solution that is compatible with Biocare's series of RISH™ probes for rapid *in situ* hybridization. The need for multiple retrieval buffers including EDTA, citrate buffer or high pH Tris buffers is eliminated when the use of RISH™ Retrieval is employed. When used in combination with RISHzyme™ for *in situ* hybridization, a synergistic effect on probe accessibility to nucleic acid targets is achieved. RISH™ Retrieval incorporates Assure™ technology, a color-coded high temperature pH indicator solution. The end-user is assured by visual inspection that the solution is at the correct dilution and pH. This product is specially formulated for superior pH stability at high temperature. RISH™ Retrieval is odorless, non-toxic, non-flammable, and sodium azide and thimerosal free.

**Known Applications:**

*in situ* hybridization (FFPE)

**Supplied As:**

**500ml**

RISH™ Retrieval, 10X concentrate (BRR0209M)

**Materials and Reagents Needed But Not Provided:**

Materials and Reagents (Needed But Not Provided):

Microscope slides, positively charged

Desert Chamber\* (Drying oven)

Positive and negative tissue controls

Xylene (Could be substituted with xylene substitute\*)

Ethanol or reagent alcohol

Decloaking Chamber\* (Pressure cooker)

Deionized or distilled water

Wash buffer\*(TBS/PBS)

RISH™ Probes\*

RISH™ Detection Kits\*

Hematoxylin\*

Bluing reagent\*

Mounting medium\*

\* Biocare Medical Products: Refer Biocare Medical Catalog for further information regarding catalog number and ordering information. Certain reagents listed above are based on specific application and detection system used.

**Storage and Stability:**

Store at room temperature. Do not use after expiration date printed on vial. If reagents are stored under conditions other than those specified in the package insert, they must be verified by the user. Diluted reagents should be used promptly; any remaining reagent should be stored at room temperature.

**Staining Protocol Recommendations:**

1. Deparaffinize tissues and hydrate to water. If necessary, block for endogenous peroxidase and wash in DI water.
2. Perform digestion according to RISH procedure.
3. Dilute concentrated RISH™ Retrieval at a ratio of 1:10 (1 ml of RISH™ to 9 ml of deionized water).
4. Place slides into 1X retrieval solution in a slide container (e.g. Coplin Jar, Tissue-Tek™ staining dish, or metal slide canister).
5. Retrieve sections using Biocare's Decloaking Chamber at 80-90°C for 15 minutes. Follow the recommendations on the RISH™ Detection kit data sheet.
6. Check solution for appropriate color change. (See Technical Note #1)

**Staining Protocol Recommendations Cont'd:**

7. Gently rinse by gradually adding DI water to the solution, then remove slides and rinse with DI water.

**Technical Notes:**

1. Concentrated RISH™ Retrieval is a bright yellow color. RTU or 1X solution is a pale yellow color. When the solution reaches 80-125°C, the solution turns yellow and indicates that the high temperature solution is at correct pH. Should the pH rise above 7.0, the solution turns a fuschia red color. Should the pH drop too low, the solution turns a pink color.
2. If using Biocare's Desert Chamber Pro (a programmable turbo-action drying oven), dry sections at 25°C overnight or at 37°C for 30-60 minutes and then dry slides at 60°C for 30 minutes.
3. Use positive charged slides (use Biocare's Kling-On HIER Slides) and cut tissues at 4-5 microns. Do not use any adhesives in the water bath. Poor fixation and processing of tissues will cause tissue sections to fall off the slides, especially fatty tissues such as breast. Tissues should be fixed a minimum of 6-12 hours.
4. Protocol time and temperatures for HIER can vary depending on the Decloaking Chamber model used. Please refer to the relevant Decloaking Chamber manual for appropriate times and temperatures.

**Limitations:**

This product is provided for Research Use Only (RUO) and is not for use in diagnostic procedures. Suitability for specific applications may vary and it is the responsibility of the end user to determine the appropriate application for its use.

**Precautions:**

1. This product is not classified as hazardous. The preservative used in this reagent is Proclin 300 and the concentration is less than 0.25%. Overexposure to Proclin 300 can cause skin and eye irritation and irritation to mucous membranes and upper respiratory tract. The concentration of Proclin 300 in this product does not meet the OSHA criteria for a hazardous substance. Wear disposable gloves when handling reagents.
2. Specimens, before and after fixation and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.
3. Microbial contamination of reagents may result in an increase in nonspecific staining.
4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any such change.
5. Do not use reagent after the expiration date printed on the vial.
6. The SDS is available upon request and is located at <http://biocare.net/>.
7. Consult OSHA, federal, state or local regulations for disposal of any toxic substances. Proclin™ is a trademark of Rohm and Haas Company, or of its subsidiaries or affiliates.

**Technical Support:**

Contact Biocare's Technical Support at 1-800-542-2002 for questions regarding this product.