



11-154-C025

## Monoclonal Antibody to CD70 Purified Antibody (0.025 mg)

<b>Clone:</b>	Ki-24
<b>Isotype:</b>	Mouse IgG3
<b>Specificity:</b>	The mouse monoclonal antibody Ki-24 recognizes CD70, an approximately 50 kDa type II transmembrane glycoprotein expressed on activated lymphocytes and some B cell leukemias. Workshop Code: IV A109
<b>Regulatory Status:</b>	RUO
<b>Species Reactivity:</b>	Human
<b>Application:</b>	Flow Cytometry
<b>Purity:</b>	> 95% (by SDS-PAGE)
<b>Purification:</b>	Purified by protein-A affinity chromatography
<b>Concentration:</b>	1 mg/ml
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
<b>Storage / Stability:</b>	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.
<b>Expiration:</b>	See vial label
<b>Lot Number:</b>	See vial label
<b>Background:</b>	CD70, also known as TNFSF7 or CD27L, is a 50 kDa type II transmembrane glycoprotein of the TNF superfamily. It is expressed mainly on activated lymphocytes, including NK cells, and forms trimeric structure. CD70 plays a role in T-cell activation, proliferation and differentiation, in enhancing the generation of cytolytic T cells, and in long-term maintenance of T cell memory. It is also involved in B cell differentiation induced by activated plasmacytoid dendritic cells, which also express CD70.
<b>References:</b>	*Kaufmann Y, Amariglio N, Rosenthal E, Hirsch YJ, Many A, Rechavi G: Proliferation response of leukemic cells to CD70 ligation oscillates with recurrent remission and relapse in a low-grade lymphoma. <i>J Immunol.</i> 2005 Nov 15;175(10):6940-7. *Trebing J, El-Mesery M, Schäfer V, Weisenberger D, Siegmund D, Silence K, Wajant H. CD70-restricted specific activation of TRAILR1 or TRAILR2 using scFv-targeted TRAIL mutants. <i>Cell Death Dis.</i> 2014 Jan 30;5:e1035. *Israel BF, Pickles RJ, Segal DM, Gerard RD, Kenney SC: Enhancement of adenovirus vector entry into CD70-positive B-cell Lines by using a bispecific CD70-adenovirus fiber antibody. <i>J Virol.</i> 2001 Jun;75(11):5215-21. Kaufmann Y, Amariglio N, Rosenthal E, Hirsch YJ, Many A, Rechavi G: Proliferation response of leukemic cells to CD70 ligation oscillates with recurrent remission and relapse in a low-grade lymphoma. <i>J Immunol.</i> 2005 Nov 15;175(10):6940-7. *Hashimoto-Okada M, Kitawaki T, Kadowaki N, Iwata S, Morimoto C, Hori T, Uchiyama T: The CD70-CD27 interaction during the stimulation with dendritic cells promotes naive CD4(+) T cells to develop into T cells producing a broad array of immunostimulatory cytokines in humans. <i>Int Immunol.</i> 2009 Aug;21(8):891-904.

**For laboratory research only, not for drug, diagnostic or other use.**



**Antibodies**

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic or therapeutic use. Not for resale or transfer either as a stand-alone product or as a component of another product without written consent of EXBIO. EXBIO will not be held responsible for patent infringement or other violations that may occur with the use of our products. All orders are accepted subject to EXBIO's term and conditions which are available at [www.exbio.cz](http://www.exbio.cz).

**For laboratory research only, not for drug, diagnostic or other use.**

---

EXBIO Praha | Nad Safinou II 341 | 252 50 Vestec u Prahy | Czech Republic  
Tel: +420 261 090 666 | Fax: +420 261 090 660 | [orders@exbio.cz](mailto:orders@exbio.cz) | [www.exbio.cz](http://www.exbio.cz)