



11-581-C100

## Monoclonal Antibody to CD326 / EpCAM Purified Antibody (0.1 mg)

<b>Clone:</b>	VU-1D9
<b>Isotype:</b>	Mouse IgG1
<b>Specificity:</b>	The mouse monoclonal antibody VU-1D9 recognizes an epitope within EGF-like domain I of CD326 / EpCAM, a marker of epithelial lineages. This antibody strongly stains various normal epithelial cells and carcinomas.
<b>Regulatory Status:</b>	RUO
<b>Immunogen:</b>	Small cell lung carcinoma cell line H69.
<b>Species Reactivity:</b>	Human, Other species Not tested
<b>Application:</b>	Flow Cytometry Recommended dilution: 1-3 µg/ml Immunohistochemistry (paraffin sections) Recommended dilution: 5-10 µg/ml Positive control: colon epithelium Application note: antigen retrieval: heat (sodium citrate) + trypsin; mAb incubation: 1 hour / RT; detection DAB Immunoprecipitation Recommended dilution: 1-2 µg / 100-500 µg of protein in 1 ml lysate Western Blotting Recommended dilution: 1-4 µg/ml Immunohistochemistry (frozen sections) Immunocytochemistry Recommended dilution: 4 µg/ml
<b>Purity:</b>	> 95% (by SDS-PAGE)
<b>Purification:</b>	Purified by protein-A affinity chromatography
<b>Concentration:</b>	1 mg/ml
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
<b>Storage / Stability:</b>	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.
<b>Expiration:</b>	See vial label
<b>Lot Number:</b>	See vial label
<b>Background:</b>	CD326 / EpCAM (also known as ESA, EGP40, EGP-2, KSA1/4, CO17-1A, GA733-2, MOC31, Ber-EP4) is a 40 kDa transmembrane glycoprotein serving as adhesion molecule in the basolateral membranes in a variety of epithelial cells. CD326 mediates calcium-independent homotypic cell-cell adhesions. CD326 over-expression has been detected in many epithelial tumours and is often associated with bad prognosis. It has been used for diagnostics of (pre-) malignancies at early stages.

**For laboratory research only, not for drug, diagnostic or other use.**



**Antibodies**

- References:**
- \*Tsubura A, Senzaki H, Sasaki M, Hilgers J, Morii S: Immunohistochemical demonstration of breast-derived and/or carcinoma-associated glycoproteins in normal skin appendages and their tumors. *J Cutan Pathol.* 1992 Feb;19(1):73-9.
  - \*Ogura E, Senzaki H, Yoshizawa K, Hioki K, Tsubura A: Immunohistochemical localization of epithelial glycoprotein EGP-2 and carcinoembryonic antigen in normal colonic mucosa and colorectal tumors. *Anticancer Res.* 1998 Sep-Oct;18(5B):3669-75.
  - \*Li G, Passebosc-Faure K, Lambert C, Gentil-Perret A, Blanc F, Oosterwijk E, Mosnier JF, Genin C, Tostain J: Flow cytometric analysis of antigen expression in malignant and normal renal cells. *Anticancer Res.* 2000 Jul-Aug;20(4):2773-8.
  - \*Winter MJ, Nagtegaal ID, van Krieken JH, Litvinov SV: The epithelial cell adhesion molecule (Ep-CAM) as a morphoregulatory molecule is a tool in surgical pathology. *Am J Pathol.* 2003 Dec;163(6):2139-48.
  - \*Brunner A, Prelog M, Verdorfer I, Tzankov A, Mikuz G, Ensinger C: EpCAM is predominantly expressed in high grade and advanced stage urothelial carcinoma of the bladder. *J Clin Pathol.* 2008 Mar;61(3):307-10.

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic or therapeutic use. Not for resale or transfer either as a stand-alone product or as a component of another product without written consent of EXBIO. EXBIO will not be held responsible for patent infringement or other violations that may occur with the use of our products. All orders are accepted subject to EXBIO's term and conditions which are available at [www.exbio.cz](http://www.exbio.cz).

**For laboratory research only, not for drug, diagnostic or other use.**

---

EXBIO Praha | Nad Safinou II 341 | 252 50 Vestec u Prahy | Czech Republic  
Tel: +420 261 090 666 | Fax: +420 261 090 660 | [orders@exbio.cz](mailto:orders@exbio.cz) | [www.exbio.cz](http://www.exbio.cz)