

1A-108-C025

## **Monoclonal Antibody to Cytokeratin (Pan-reactive)** Allophycocyanin (APC) conjugated (0.025 mg)

Clone: C-11

Mouse IgG1 Isotype:

Specificity: The antibody C-11 reacts with Cytokeratin peptides 4, 5, 6, 8, 10, 13, 18.

Cytokeratins are a member of intermediate filaments subfamily represented in

epithelial tissues.

**Regulatory Status: RUO** 

Immunogen: Keratin-enriched preparation from human epidermoid carcinoma cell line A431.

**Species Reactivity:** Mammalian

Preparation: The purified antibody is conjugated with cross-linked Allophycocyanin (APC) under

optimum conditions. The conjugate is purified by size-exclusion chromatography.

0.1 mg/ml **Concentration:** 

The reagent is provided in stabilizing phosphate buffered saline (PBS) solution Storage Buffer:

containing 15mM sodium azide.

Storage / Stability: Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not

use after expiration date stamped on vial label.

The reagent is designed for Flow Cytometry analysis. Usage:

> Suggested working dilution is 5 µg/ml. Indicated dilution is recommended starting point for use of this product. Working concentrations should be determined by the

investigator.

See vial label **Expiration:** 

See vial label Lot Number:

Cytokeratins are a subfamily of intermediate filaments and characterized by **Background:** 

remarkable biochemical diversity. Cytokeratins are represented in epithelial tissues by at least 20 different polypeptides, molecular weight between 40 kDa and 68 kDa. The individual cytokeratin polypeptides are designated 1 to 20 and divided into the type I (acidic cytokeratins 9-20) and type II (basic to neutral cytokeratins

1-8) families.

\*Kovarik J, Rejthar A, Lauerova L, Vojtesek B, Bartkova J: Monoclonal antibodies References:

against individual cytokeratins in the detection of metastatic spread. Int J Cancer

Suppl. 1988;3:50-5.

\*Vojtĕsek B, Stasková Z, Nenutil R, Lauerová L, Kovarík J, Rejthar A, Bártková J, Bártek J: Monoclonal antibodies recognizing different epitopes of

cytokeratin No.18. Folia Biol (Praha). 1989;35(6):373-82.

\*Bartek J, Vojtesek B, Staskova Z, Bartkova J, Kerekes Z, Rejthar A, Kovarik J: A series of 14 new monoclonal antibodies to keratins: characterization and value in

diagnostic histopathology. J Pathol. 1991 Jul;164(3):215-24. \*Hamakawa H, Sumida T, Tanioka H, Sogawa K, Yamada T: Extraction of cytokeratin from the human submandibular gland and its electrophoretic analysis.

Res Commun Mol Pathol Pharmacol. 1998 Aug;101(2):115-26.

\*Broekema M, Harmsen MC, Koerts JA, Petersen AH, van Luyn MJ, Navis G, Popa ER: Determinants of tubular bone marrow-derived cell engraftment after

renal ischemia/reperfusion in rats. Kidney Int. 2005 Dec;68(6):2572-81.

For laboratory research only, not for drug, diagnostic or other use.



## PRODUCT DATA SHEET

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic or therapeutic use. Not for resale or transfer either as a stand-alone product or as a component of another product without written consent of EXBIO. EXBIO will not be held responsible for patent infringement or other violations that may occur with the use of our products. All orders are accepted subject to EXBIO's term and conditions which are available at www.exbio.cz.

For laboratory research only, not for drug, diagnostic or other use.