

Data Sheet

Strep-tag II peptide

Cat. No.: 2-1018-002; 2-1018-009

Lot No.: 1018 -

Last date of revision
March 2012
Version 1018-3

Description	Synthetic <i>Strep-tag II</i> peptide recognized by <i>Strep-Tactin</i>[®] and the monoclonal <i>Strep-tag</i>[®] antibodies <i>StrepMAB-Classic</i> and <i>StrepMAB-Immo</i> Sequence: SAWSHPQFEK, MW=1216
Application	Can be used for competitive elution of <i>Strep-tag II</i> proteins bound to <i>StrepMAB-Classic</i> columns; Cat. No. 2-1526-001 or Cat. No. 2-1526-501). Working concentration: 0.5 mM Not suitable for elution of proteins immobilized on <i>StrepMAB-Immo</i> .
Use	Dissolve the whole amount in 3 ml suitable buffer to get a 0.5 mM solution. It is recommended to use the same buffer as used for washing the column.
Form	Lyophilized as TFA salt.
Amount	1.8 mg peptide (Cat. No. 2-1018-002) 5 x 1.8 mg peptide (Cat. No. 2-1018-009)
Stability	6 months after shipping
Storage	-20°C
Shipment	Room temperature

For research use only

Important licensing information

This product is based on *Strep-tag* technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.