

Please note this data sheet has been changed effective October 24, 2014

Phosphate Buffered Saline

INTENDED USE

For In Vitro Diagnostic Use

<u>AVAILABILITY:</u>	<u>Catalog #</u>	<u>Volume</u>
	AP-9009-10	powder for 1L of 10X

SPECIFICITY: N/A

ENZYME: N/A

CHROMOGEN/SUBSTRATE: N/A

DESCRIPTION

The packet contains a 10X concentrated formulation of phosphate buffered saline that can be used as a wash buffer during immuno-labeling procedures.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature. When stored at room temperature, it is stable for three years.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

To make a 10X concentrated phosphate buffered saline empty the contents of the pack into a container and add distilled deionized water to the final volume of 1lit. This is a 10X stock solution and should be diluted 10-fold with distilled water to make the 1X working solution.

REFERENCES

N/A

TROUBLESHOOTING

Please contact Thermo Fisher Scientific Technical Support by phone (1-269-544-5600 or 1-800-522-7270) or by email (lab.reagents@thermofisher.com).

