

Product Datasheet

Bladder Lysate NBP2-29005

Unit Size: 4 Vials

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



support@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-29005

Updated 1/9/2014 v.20.1

NBP2-29005

Human Bladder Tissue Lysate (Adult)

| Product Information | |
|------------------------------------|--|
| Unit Size | 4 Vials |
| Concentration | 1 mg/ml |
| Storage | Store at -80C. Avoid freeze-thaw cycles. |
| Reconstitution Instructions | The vial is provided with a 10% overfill. Maximum recovery can be obtained by centrifuging the vial briefly to collect any solution on the cap and tube sides. |
| Product Description | |
| Species | Human |
| Specificity/Sensitivity | Age :57 Diagnosis:Transitional cell carcinoma Stage:I TNM:T1bN0Mx |
| Immunogen | Clinical Tissue Human Bladder Tissue Lysate Bundle Set |
| Preparation Method | Tissue specimens are homogenized in modified RIPA buffer to obtain the soluble proteins, and centrifuged to clarify. The pellet was further extracted with a second buffer to obtain the less soluble protein fraction. The lysate solution may appear turbid at cold temperatures due to insolubility of buffer components. The solution should clear upon warming to room temperature. |
| Lysate Type | Tissue |
| Lysate Tissue | Bladder |
| Lysate Life Stage | Adult |
| Product Application Details | |
| Application Notes | These lysates are proteomic discovery tools.Researchers should validate and optimize for individual use. Potential applications MAY include WB, immunoprecipitation, protein-protein interactions, ligand binding, ELISA. Note: For use in 1D and 2D gel electrophoresis, the addition of a denaturing gel loading buffer with reducing agents may be required. |



Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP2-29005

| | |
|------------|---|
| NBP2-30252 | Bladder Tissue Slides (Transitional Cell Carcinoma) |
|------------|---|

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our guarantee, please visit www.novusbio.com/guarantee.

