

Product datasheet for CV500001

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

AAV Serotyping Kit I

Product data:

Product Type: AAV Control Particles

Description: This kit contains pre-made AAV-2, AAV-5, AAV-6, AAV-8, and AAV-9 control particles for testing

each serotype's infection efficiency in your desired cell type.

Promoter: CMV

Validation: To validate the quality of our AAV control particles, we have developed the following

procedure:

AAV genome titer is determined by qPCR. AAV purity is determined by Silver staining.

AAV transduction efficiency is analyzed by in vitro transduction in HEK293T.

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at

-20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during

the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible,

particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles.