

BioMag Solo-Sep Microcentrifuge Tube Separator

Description

The Solo-Sep is a durable, plastic magnetic separation unit designed to accommodate one 1.5ml microcentrifuge tube. The magnet and tube holder are positioned to efficiently magnetically separate microliter volumes common in molecular biology applications. This unit will separate most 40 ul to 1ml volumes containing <5mg of magnetic particles.

Use

Magnetic separations can be performed on one 1.5ml microcentrifuge tube. Once the tube is placed in the tube holder, the magnetic particles are pulled to the side of the unit to complete a separation in less than one minute. Upon completion of the separation, the supernatant may be aspirated, leaving the magnetic particles as a pellet on the side of the tube.

Care

The BioMag Solo-Sep unit should be wiped clean with mild detergent.

DO NOT IMMERSE. Harsh chemicals and organic solvents may cause damage to the surface of the unit.

Warning

Strong Magnetic Field. This separation unit should be used and stored at a safe distance from magnetic media such as computers, videotapes and credit cards etc.

Ordering Information:

Cat. #	Description	Size
8MB4112S	BioMag Solo-Sep Microcentrifuge Tube Separator	1 ea

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

