

Salusin-alpha (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-010-67	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

Salusin

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-010-67

Standard Size 96 well


Sensitivity Minimum Detectable Concentration = 44.3 pg/ml

Range 0-10000 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Salusin-alpha (H)	100
	Salusin-alpha-Gly29 (H)	100
	Salusin (50-mer) (H)	100
	Salusin-beta (H)	0
	Salusin (40-mer) (H)	0
	Endothelin-1 (H, R)	0

Sequence Ser - Gly - Ala - Leu - Pro - Pro - Ala - Pro - Ala - Ala - Pro - Arg - Pro - Ala - Leu - Arg - Ala - Gln - Arg - Ala - Gly - Pro - Ala - Gly - Pro - Gly - Ala - Lys - NH₂

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 47.4-1181 pg/ml

ED50 212 pg/ml

Sample Level **65 pg/ml**