

Magnetic Luminex® Performance Assay Human IL-1β/IL-1F2 Kit

Catalog Number: LUHM201
Pack Size: 100 Tests

SPECIFICATIONS AND USE

Recommended Sample Types Microparticle Region Components

- Cell culture supernates, serum, EDTA plasma, and heparin plasma.
- Region-20
- Microparticle Concentrate (Part 894430) is supplied as a 100X concentrated stock (0.075 mL) with preservatives.
- Biotin-Antibody Concentrate (Part 892618) is supplied as a 100X concentrated stock solution (0.075 mL) with preservatives.

Other Supplies Required

Storage


- Magnetic Luminex Performance Assay Human Base Kit A (Catalog Number LUHM000).
- Store the unopened kit at 2-8 °C. Do not use past the expiration date on the label.
- Avoid freezing microparticles.
- Protect microparticles from light.

Instructions for Use

• Refer to the Base Kit insert for the Luminex Performance Assay procedure.

TYPICAL DATA

This human IL-1 β standard curve is provided only for demonstration. A standard curve must be generated each time an assay is run, utilizing values from the Standard Value Card included in the Base Kit.



Standard	pg/mL	MFI	Average	Corrected
Blank	0	12 13	13	
1	2400	10,168 10,272	10,220	10,207
2	800	4662 4684	4673	4660
3	267	1858 1889	1874	1861
4	89	817 834	825	812
5	30	345 351	348	335
6	10	144 148	146	133
7	3	72 73	73	60

PERFORMANCE CHARACTERISTICS

All data were collected with assays run as a multiplex.

Data obtained with polystyrene and magnetic beads were equivalent.

Sensitivity - The Minimum Detectable Dose (MDD) was determined by adding two standard deviations to the MFI of twenty zero standard replicates and calculating the corresponding concentration.

Forty-two assays were evaluated, and the MDD of human IL-1β ranged from 0.11-0.57 pg/mL. The mean MDD was 0.27 pg/mL.

752771.0

Intra-assay Precision (precision within an assay) - Three samples of known concentration were tested twenty times on one plate to assess precision within an assay.

Inter-assay Precision (precision between assays) - Three samples of known concentration were tested in twenty separate assays to assess precision between assays.

	Int	Intra-assay Precision				Inter-assay Precision		
Sample	1	2	3		1	2	3	
n	20	20	20		20	20	20	
Mean (pg/mL)	29	95	662		34	115	713	
Standard Deviation	1.6	4.3	33.5		3.4	11.3	52.5	
% CV	5.5	4.5	5.1		10.0	9.8	7.4	

Recovery and Linearity – Samples containing and/or spiked with high concentrations of IL-1 β were evaluated for recovery and were serially diluted to evaluate assay linearity.

Recovery					
Sample Type	Average % Recovery	Range (%)			
Cell culture supernates	93	92-94			
Serum	109	98-119			
EDTA plasma	100	88-111			
Heparin plasma	99	87-106			

Linearity							
		Cell culture supernates	Serum	EDTA Plasma	Heparin Plasma		
1:2	Average % of Expected	99	103	105	99		
1.2	Range (%)	89-116	99-111	87-115	92-104		
1.4	Average % of Expected	94	101	106	98		
1:4	Range (%)	79-112	95-104	91-114	90-104		
1:8	Average % of Expected	99	100	113	97		
	Range (%)	75-119	92-110	103-121	88-101		
			•		•		

Specificity - This assay recognizes natural and recombinant human IL-1 β . The assay was tested for cross-reactivity and interference with the following factors. Less than 0.5% cross-reactivity and interference was observed.

Recombinant human:			Recombinant mouse:		Recombinant rat:	Recombinant porcine:	Recombinant human multiplex partners:	
6Ckine	IL-1 RII	IL-17	G-CSF	IL-8	GM-CSF	GM-CSF	ENA-78	IL-8
CNTF	IL-2 R $lpha$	IL-18	GM-CSF	IL-10	IFN-γ	IL-1 α	FGF basic	IL-10
β-ECGF	IL-2 Rβ	LIF	IFN-γ	IL-17	IL-1α	IL-1β	G-CSF	IL-17
FGF acidic	IL-2 Rγ	LIF R	IL-1α	MIP-1 $lpha$	IL-1β	IL-2	GM-CSF	MCP-1
FGF-4	IL-3 Rα	MIP-1 α	IL-1ra	MIP-1β	IL-2	IL-4	IFN-γ	MIP-1 α
FGF-5	IL-4 R	MIP-3 α	IL-1	RANTES	IL-4	IL-6	IL-1α	MIP-1β
FGF-6	IL-5 R $lpha$	MIP-3β	IL-2	Тро	IL-6	IL-8	IL-1ra	RANTES
FGF-9	IL-6 R	MCP-2	IL-4	TNF-α	IL-10	IL-10	IL-2	Тро
FGF-10	IL-10 R	MCP-3	IL-5	VEGF	TNF- α	Leptin	IL-4	TNF-α
FGF-18	IL-3	MCP-4	IL-6			TNF-α	IL-5	VEGF
GCP-2	IL-7	M-CSF					IL-6	
$GRO\alpha$	IL-9	TNF RI						
GR0β	IL-11	TNF- $lpha$						
GROγ	IL-12 p40	VEGF ₁₂₁						
I-309	IL-12 p70	VEGF ₁₆₅						
IGF-I	IL-13	VEGF-D						
IGF-II	IL-15							
IL-1 RI	IL-16							

TECHNICAL HINTS

- Protect the microparticles and streptavidin-PE from light at all times.
- Refer to the Base Kit Standard Value Card for reconstitution volume and values of the reconstituted standard.
- Diluted microparticles cannot be stored. Make a fresh dilution of microparticles each time the assay is run.
- The use of a magnetic device made to accommodate a microplate is necessary for washing.
- Discrepancies may exist in values obtained for the same analyte utilizing different technologies.

Luminex Performance Assays afford the user the benefit of multianalyte analysis of biomarkers in a complex sample. For each sample type, a single, multipurpose diluent is used to optimize recovery, linearity, and reproducibility. Such a multipurpose diluent may not optimize any single analyte to the same degree that a unique diluent selected for analysis of that analyte can optimize conditions. Therefore, some performance characteristics may be more variable than those for assays designed specifically for single analyte analysis.

752771.0 www.RnDSystems.com