

# TREVIGEN<sup>®</sup> Product Data

For Research Use Only. Not For Use In Diagnostic Procedures.

---

## TACS•XL<sup>™</sup>-Replenisher (apoTACS-replenisher)

Catalog #: 4828-30-R

**Description:** Replenisher reagents for TACS•XL-Basic (Cat# 4828-30-K), TACS•XL-Blue Label (Cat# 4828-30-BK) or TACS•XL-DAB (Cat# 4828-30-DK).

### Contents:

<u>Component</u>	<u>Quantity</u>	<u>Storage</u>	<u>Catalog #</u>
Proteinase K Solution	30 µl	-20°C	4800-30-01
B-dNTP Mix	30 µl	-20°C	4828-30-04
TdT Enzyme	30 µl	-20°C	4810-30-05
Strep-HRP	30 µl	4°C	4800-30-06
Anti-BrdU Antibody	30 µl	4°C	4828-30-06
Strep-Diluent	7.5 ml	4°C	4828-30-12

**Storage:** Caution! There are multiple storage temperatures for these products. Please follow appropriate storage temperatures as listed above and on the individual product labels.

**Applications:** Replenisher reagents for the *in situ* detection of apoptosis using TACS•XL.

**Note:** For your convenience, we have included the TACS•XL-DAB as well as the TACS•XL-Blue Label complete protocols with your TACS•XL-Replenisher reagents. Please select the appropriate procedure before starting your experiment. If you have any questions, please call our technical service department.

© 2001 Trevigen, Inc. All Rights Reserved. Trevigen is a registered trademark and TACS•XL and Cytonin are trademarks of Trevigen, Inc. v01121

---

## TREVIGEN<sup>®</sup>

8405 Helgerman Court, Gaithersburg, MD 20877 USA  
Voice: 1-800-TREVIGEN (1-800-873-8443) • 301-216-2800  
Fax: 301-216-2801 • e-mail: [info@trevigen.com](mailto:info@trevigen.com) • [www.trevigen.com](http://www.trevigen.com)