



---

## CERTIFICATE OF ANALYSIS

Product **Sialic Acid Agarose**

Catalog No. AC-1008

Amount 10 ml of settled gel

Lot No. ZA0605

Matrix: Cross-linked, 4%, 90 $\mu$ -120 $\mu$  diameter agarose beads

Linkage: Randomly coupled through a proprietary 5-atom hydrophilic linker arm

Suspension Solution HBS, pH 7.5, 0.08% azide

Binding Capacity Binds 5-6 mg *Sambucus nigra* per ml of gel

Storage Conditions Refrigerate - DO NOT FREEZE

### Typical Protocol for Binding and Elution

- 1) Wash with two column volumes of HBS, pH 7.5, or TBS, pH 7.5
- 2) Pass the solution to be fractionated over the column. (Note: if binding is poor at room temperature, equilibrate column and extract at 4 °C prior to the column pass.)
- 3) Wash the column with the equilibrating buffer at the temperature at which the pass was made until the OD<sub>280nm</sub> absorbance of the wash decreases to background.
- 4) Elute the bound protein with Eluting Solution (Cat. No. ES-7100).