

## Visikol® 1.75mm ClearWell Silicon Imaging Chambers

### Microscope Slide Imaging Chamber

Catalog Number – CW1.75

#### INTRODUCTION

After you have finished staining and clearing your tissues, you can use ClearWells to mount your tissue for imaging. ClearWells can be used with any tissue clearing technique, are chemically inert, and naturally adhere to glass microscope slides. Simply place a cover slip on top and start imaging.


#### TYPICAL APPLICATIONS

- Whole mount microscope slide sample prep
- Archiving of samples
- Compatible with confocal, multiphoton, widefield, and traditional light microscopy

#### ADVANTAGES

- Compatible with all Visikol HISTO reagents
- Creates mounted slides which can be stored long term

#### PRODUCT IMAGE



#### DISCLAIMER

No freedom from infringement of any patent owned by Visikol or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer's use and for ensuring that the Customer's workplace and disposal practices are in compliance with applicable laws and other governmental enactments. The product shown in this literature may not be available for sale and/or available in all geographies where Visikol is represented. The claims made may not have been approved for use in all countries. Visikol assumes no obligation or liability for the information in this document. **NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.**