



## Visikol® HISTO-1™

### Optical Clearing Agent

#### INTRODUCTION

Visikol HISTO-1 is an optical clearing agent for biological tissues enabling rapid and easy visualization of tissues in 3D without altering the underlying cellular morphology. Visikol HISTO-1 is fully compatible with immunostaining as well as many common other common fluorescent labels. Using Visikol HISTO-1 with confocal or light sheet microscopy allows researchers to see deeper into tissues than ever before.

#### TYPICAL APPLICATIONS

- Optical clearing agent in tissue imaging studies
- Refractive index matching agent for microscopy

#### ADVANTAGES

- Easy-to-use tissue clearing modality
- Rapid tissue clearing
- Clearing is completely reversible with no alteration to tissue

#### TYPICAL PHYSICAL PROPERTIES

These properties are typical but should not be considered specifications. Please contact us to obtain the certificate of analysis for your specific lot of material.

Appearance:	Clear, colorless liquid
Refractive Index:	1.49-1.50
Water:	0.05-0.25%
Maximum Absorbance 450-750nm:	0.01-0.10 AU
Maximum Fluorescence Emission 450-750nm:	1-10 IU

#### DISCLAIMER

No freedom from infringement of any patent owned by Visikol or others is to be inferred. Because use conditions and applicable laws may differ from one location to another and may change with time, the Customer is responsible for determining whether products and the information in this document are appropriate for the Customer’s use and for ensuring that the Customer’s workplace and disposal practices are in compliance with applicable laws and other governmental enactments. The product shown in this literature may not be available for sale and/or available in all geographies where Visikol is represented. The claims made may not have been approved for use in all countries. Visikol assumes no obligation or liability for the information in this document. NO WARRANTIES ARE GIVEN; ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE EXPRESSLY EXCLUDED.