

Stem Cells Generation

5-Azacytidine

DNA hypomethylation reagent that increases stem cell reprogramming efficiency.¹
Enhances the efficiency of iPS cell generation.²

Product No: 10-1267

50 mg

250 mg

Bay K8644

Enables reprogramming of Act3/Klf4-transduced embryonic fibroblasts in combination with BIX-01294.³

Product No: 10-4518

5 mg

25 mg

Cordycepin

AMPK activator that maintains stem cell pluripotency and increases iPS cell generation efficiency.⁴

Product No: 10-2408

10 mg

50 mg

DZNep

Methyltransferase inhibitor that can be used in chemical cocktails to induce iPSs from somatic cells.⁵

Product No: 10-1524

1 mg

5 mg

Na Valproate

Improves stem cell reprogramming efficiency and enables efficient induction of pluripotency without introduction of the oncogene c-Myc.⁶

Product No: 10-1009

5 g

PD 325901

Induces generation of ground state iPS cells from human and rat somatic cells in combination with CHIR-99021 and A83-01.⁷

Product No: 10-2161

5 mg

25 mg

Pifithrin α

p53 inhibitor that induced pluripotency in primary bone marrow derived mesenchymal stromal cells in conjunction with SB431542 and PD325901.⁸

Product No: 10-2480

5 mg

25 mg

Purmorphamine

Part of a chemical cocktail that reprograms mouse somatic cells into pluripotent stem-like cells⁹ and human astroglial cells into functional neurons¹⁰.

Product No: 10-4622

5 mg

25 mg

RepSox

Potent ALK5 inhibitor that can replace Sox2 in reprogramming adult cells into pluripotent stem cells.¹¹

Sulindac Sulfone No: 10-4546

5 mg

25 mg

Thiazovivin

Dramatically improves the efficiency of induced pluripotent stem cell generation from human fibroblasts.¹²

Product No: 10-1191

5 mg

25 mg


Tranylcypromine hemisulfate

Induces reprogramming of human primary keratinocytes into iPS cells in conjunction with CHIR99021.¹³


Product No: 10-4535

50 mg

250 mg



Purmorphamine



Thiazovivin

REFERENCES

1. Mikkelsen *et al.* (2008) *Nature* **454** 49
2. De Carvalho *et al.* (2010) *Trends Cell Biol.* **20** 609
3. Shi *et al.* (2008) *Cell Stem Cell* **3** 568
4. Wang *et al.* (2020) *Sci. Rep.* **10** 2187
5. Hou *et al.* (2013) *Science* **341** 651
6. Hangfu *et al.* (2008) *Nat. Biotechnol.* **26** 795
7. Li *et al.* (2009) *Cell Stem Cell* **4** 16
8. Megges *et al.* (2015) *Stem Cell Res.* **15** 565
9. Kang *et al.* (2015) *Biomaterials* **35** 7336
10. Zhang *et al.* (2015) *Cell Stem Cell* **17** 735
11. Ichida *et al.* (2009) *Cell Stem Cell* **5** 491
12. Lin *et al.* (2009) *Nat. Methods* **6** 805
13. Li *et al.* (2009) *Stem Cells* **27** 2992
14. Wang *et al.* (2019) *Cell Biol.* **97** 702