

# Atlas of PROGEN's Keratin Antibodies



**Lung:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)



**Prostate:** K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)



**Kidney:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)



**Colon:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429), K20 (65054, 61033, GP-K20, 61032, 65026, 65133, 61054, 65132, 65126, 65154, 65133, 61454)




**Ovary:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)

**Cervix:** K4 (GP-K4), K5 (GP-CK5), K6 (GP91, 61090, GP-K6, 65190), K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K10 (11414, GP-K10), K13 (GP-K13, 61007, 65007), K14 (GP-CK14), K15 (GP-CK15, GP-K15), K16 (GP-K16), K17 (65136, GP-CK17, 61036), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)

**Uterus:** K5 (GP-CK5), K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)

**Breast:** K5 (GP-CK5), K6 (GP91, 61090, GP-K6, 65190), K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K14 (GP-CK14), K17 (65136, GP-CK17, 61036), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)



**Mesothelioma:** K5 (GP-CK5), K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)

**Biliary Tract:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429), K20 (65054, 61033, GP-K20, 61032, 65026, 65133, 61054, 65132, 65126, 65154, 65133, 61454)

**Pancreas:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429), K20 (65054, 61033, GP-K20, 61032, 65026, 65133, 61054, 65132, 65126, 65154, 65133, 61454)

**Skin:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429), K20 (65054, 61033, GP-K20, 61032, 65026, 65133, 61054, 65132, 65126, 65154, 65133, 61454)



**Liver:** K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429)



**Stomach:** K7 (61025, 65025, GP-K7), K8 (65130, 61038, GP-K8, 65138, 61538), K18 (61028, 65134, GP-K18, 65137, 65128, GP-CK18, 61034, 65028, 11416, 61528, 61428), K19 (61010, 65129, 65010, 11417, GP-CK19, 61029, 65135, 61429), K20 (65054, 61033, GP-K20, 61032, 65026, 65133, 61054, 65132, 65126, 65154, 65133, 61454)