

Tube #	Cryoprotectant	M _r	CAS #	Formula	Class
1.	100% (+/-)-2-Methyl-2,4-pentanediol	118.18	107-41-5	CH ₃ CH(OH)CH ₂ C(CH ₃) ₂ OH	Organic, Non-volatile
2.	6.0 M 1,6-Hexanediol	118.18	629-11-8	HO(CH ₂) ₆ OH	Organic, Non-volatile
3.	100% 1,2-Propanediol	76.09	57-55-6	CH ₃ CH(OH)CH ₂ OH	Organic, Non-volatile
4.	100% 2,3-Butanediol	90.12	513-85-9	CH ₃ CH(OH)CH(OH)CH ₃	Organic, Non-volatile
5.	50% w/v NDSB-201	201.26	15471-17-7	C ₈ H ₁₁ NO ₃ S	Non-detergent
6.	6.0 M L-Proline	115.13	147-85-3	C ₅ H ₉ NO ₂	Osmolyte
7.	4.0 M Trimethylamine N-oxide dihydrate	111.14	62637-93-8	(CH ₃) ₃ NO • 2H ₂ O	Osmolyte
8.	100% Glycerol	92.09	56-81-5	HOCH ₂ CH(OH)CH ₂ O	Polyol
9.	100% Ethylene glycol	62.07	107-21-1	HOCH ₂ CH ₂ OH	Polyol
10.	50% v/v Diethylene glycol	106.12	111-46-6	(HOCH ₂ CH ₂) ₂ O	Polyol
11.	100% Polyethylene glycol 200	200	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer / Polyol
12.	100% Polyethylene glycol 400	380 - 420	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer / Polyol
13.	100% Polyethylene glycol monomethyl ether 550	550	9004-74-4	CH ₃ O(CH ₂ CH ₂ O) _n H	Polymer / Polyol
14.	80% v/v Polyethylene glycol 600	570 - 630	25322-68-2	H(OCH ₂ CH ₂) _n OH	Polymer / Polyol
15.	50% w/v Polyethylene glycol 1,000	950 - 1050	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer
16.	50% w/v Polyethylene glycol 3,350	3015 - 3685	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer
17.	50% w/v Polyethylene glycol 4,000	3,500 - 4,500	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer
18.	50% w/v Polyethylene glycol monomethyl ether 5,000	5,000	9004-74-4	CH ₃ O(CH ₂ CH ₂ O) _n H	Polymer
19.	50% w/v Polyethylene glycol 8,000	7,000 - 9,000	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer
20.	50% w/v Polyethylene glycol 10,000	8,500 - 11,500	25322-68-3	H(OCH ₂ CH ₂) _n OH	Polymer
21.	50% w/v Polyvinylpyrrolidone K 15	10,000	9003-39-8	(C ₆ H ₉ NO) _n	Polymer
22.	50% v/v Pentaerythritol propoxylate (5/4 PO/OH)	426	9051-49-4	C[CH ₂ [OCH ₂ CH(CH ₃) ₁ OH] ₄ N~ 5	Polymer
23.	100% Polypropylene glycol P 400	N/A	25322-69-4	N/A	Polymer
24.	100% Dimethyl sulfoxide (DMSO)	78.13	67-68-5	(CH ₃) ₂ SO	Solvent

Storage:

CryoPro reagents may be stored at -20 to 4°C.

Allow the reagents to return to working temperature and mix well before use.

*CryoPro contains 48 unique water soluble cryoprotectants for cryocrystallography.
Reagents formulated in Type 1+ ultrapure grade water, no pH adjustment except 39 and 44 titrated using HCl.*



Solutions for Crystal Growth

Tube #	Cryoprotectant	M _r	CAS #	Formula	Class
25.	70% w/v D-(+)-Sucrose	342.30	57-50-1	C ₁₂ H ₂₂ O ₁₁	Sugar
26.	70% w/v D-Sorbitol	182.17	50-70-4	C ₆ H ₁₄ O ₆	Sugar
27.	30% w/v Maltose monohydrate	342.30	69-79-4	C ₁₂ H ₂₂ O ₁₁ • H ₂ O	Sugar
28.	35% w/v meso-Erythritol	122.12	149-32-6	HOCH ₂ [CH(OH)] ₂ CH ₂ OH	Sugar
29.	70% w/v Xylitol	152.15	87-99-0	HOCH ₂ [CH(OH)] ₃ CH ₂ OH	Sugar
30.	15% w/v myo-Inositol	180.16	87-89-8	C ₆ H ₁₂ O ₆	Sugar
31.	20% w/v D-(+)-Raffinose pentahydrate	594.51	17629-30-0	C ₁₈ H ₃₂ O ₁₆ • 5H ₂ O	Sugar
32.	50% w/v D-(+)-Trehalose dihydrate	378.33	6138-23-4	C ₁₂ H ₂₂ O ₁₁ • 2H ₂ O	Sugar
33.	70% w/v D-(+)-Glucose monohydrate	198.17	14431-43-7	C ₆ H ₁₂ O ₆ • H ₂ O	Sugar
34.	5.0 M Lithium acetate dihydrate	102.02	6108-17-4	CH ₃ COOLi • 2H ₂ O	Salt
35.	10.0 M Lithium chloride	42.39	7447-41-8	LiCl	Salt
36.	4.0 M Lithium formate monohydrate	69.97	6108-23-2	HCO ₂ Li • H ₂ O	Salt
37.	8.0 M Lithium nitrate	68.95	7790-69-4	LiNO ₃	Salt
38.	2.0 M Lithium sulfate monohydrate	127.96	10102-25-7	Li ₂ SO ₄ • H ₂ O	Salt
39.	3.4 M Sodium malonate pH 7.0	104.06	141-82-2	C ₃ H ₄ O ₄	Salt
40.	3.4 M Magnesium acetate tetrahydrate	214.46	16674-78-5	Mg(CH ₃ COO) ₂ • 4H ₂ O	Salt
41.	5.0 M Sodium chloride	58.44	7647-14-5	NaCl	Salt
42.	7.0 M Sodium formate	68.01	141-53-7	CHNaO ₂	Salt
43.	7.0 M Sodium nitrate	84.99	7631-99-4	NaNO ₃	Salt
44.	100% Tacsimate™ pH 7.0	N/A	N/A	N/A	Salt
45.	1.0 M Sodium sulfate decahydrate	322.20	7727-73-3	Na ₂ SO ₄ • 10H ₂ O	Salt
46.	50% v/v Ethylene glycol, 25% w/v NDSB-201	N/A	N/A	N/A	Mixture
47.	30% w/v Polyethylene glycol 3,350, 20% v/v Glycerol	N/A	N/A	N/A	Mixture
48.	30% v/v Polyethylene glycol 400, 10% w/v Polyethylene glycol 20,000, 5% v/v Glycerol, 5% w/v NDSB-201	N/A	N/A	N/A	Mixture

Storage:

CryoPro reagents may be stored at -20 to 4°C.

Allow the reagents to return to working temperature and mix well before use.

*CryoPro contains 48 unique water soluble cryoprotectants for cryocrystallography.
Reagents formulated in Type 1+ ultrapure grade water, no pH adjustment except 39 and 44 titrated using HCl.*



Solutions for Crystal Growth

34 Journey
Aliso Viejo, CA 92656-3317 U.S.A.
Tel: (949) 425-1321 • Fax: (949) 425-1611
e-mail: tech@hrmail.com
Website: www.hamptonresearch.com

© 1991 - 2013 Hampton Research Corp. all rights reserved
Printed in the United States of America. This guide or
parts thereof may not be reproduced in any form without
the written permission of the publishers.