BioDynamics Laboratory Inc.

PRODUCT INFORMATION

For Research Use Only

Product Name: Sonicated Salmon Sperm DNA

Code No.: F012, F013

Size: 1 ml / tube F012 1 ml

F013 5 ml $(1ml \times 5)$

Concentration: 10 mg/ml

Description: Salmon Sperm DNA was shared by sonication.

DNA is dissolved in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Quality Control: Intrinsic DNase activity was not detected in the product,

by the use of pUC18 plasmid DNA as a control.

Preparation: Heat at 95°C for 3 minutes and cool on ice before use.

Storage Condition: - 20 °C

Related Products:

DM253	DynaMarker® Prestain Marker for Small RNA Plus
	Prestained Molecular Weight Marker for RNA (20-100 bases)
DM260	DynaMarker® Prestain Marker for RNA High
	Prestained Molecular Weight Marker for RNA (200-8,000 bases)