


TubeSpin® Bioreactor 50 - 20

DESCRIPTION

50 ml tubes with membrane filter TubeSpin® Bioreactor 50, TPP® 20 pcs.

SPECIFICATIONS

Number of tubes in set	20
Volume	50 ml



CAT. NUMBER

BS-010158-AK	50 ml tubes with membrane filter TubeSpin® Bioreactor 50, TPP® 20 pcs.
--------------	--

ACCESSORIES



RTS-1

Including TPP TubeSpin® Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1 is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-1C

Including TPP TubeSpin® Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1C is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)