



Antibody Sampler Kit  
Purified Rabbit Polyclonal Antibody (Pab) Sampler Kit

Catalog #	Applications:	Reactivity:	Accessions:
AS1009a			
Concentration:	Size:	Isotype:	Clone Name:
0.25 mg/ml	0.02 mg/antibody	Rabbit Ig	N/A
Application Data:		Calculated MW: Da Da	

Other Names:

Target/Specificity:

Select five or more antibodies from our catalog portfolio.

Application Notes:

Antibodies have been broadly tested in ELISA, Western blot, and IHC applications. Immunofluorescence data is also demonstrated for select products. Check individual product technical data sheets for details.

Format:

Purified polyclonal antibodies supplied in PBS with 0.09% (W/V) sodium azide. The antibodies are purified through protein G columns, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.

Storage:

Maintain refrigerated at 2-8 ° C for up to 6 months. For long term storage store at -20 ° C.

Precautions:

This product is for research use only. Not for use in diagnostic or therapeutic procedures.