

Goat IgG-heavy and light chain Antibody

Affinity Purified

Alkaline Phosphatase conjugated

Produced in Rabbit

Catalog No.

A50-100AP



Storage: 2 - 8 C
Amount: 1 ml at 1 mg/ml
Buffer: 50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl₂, 0.1 mM ZnCl₂ containing 0.09% NaN₃
Shelf Life: 1 year from date of receipt

Production Procedures: The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to alkaline phosphatase (alkphos).

Antibody concentration was determined by extinction coefficient prior to conjugation: absorbance at 280 nm of 1.4 equals 1.0 mg of IgG. Molar enzyme:antibody protein ratio is 1:1.

By immunoelectrophoresis and ELISA this antibody reacts specifically with goat IgG and with light chains common to other goat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species.

Country of Origin: USA

Notes: For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

This product may not be resold or modified for resale without the prior written approval of Bethyl Laboratories, Inc.

Applications: Not all listed applications have been specifically tested. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilutions immediately before use. Suggested starting dilutions are:

ELISA: 1:1000 – 1:10,000

Western Blot: 1:1000 – 1:10,000

Immunohistochemistry: 1:200 – 1:1000

Immunocytochemistry: 1:200 – 1:1000

In some cases, the antibody may be diluted further than indicated.

Related Products:	Rabbit anti-Goat IgG-heavy and light chain Antibody Affinity Purified	A50-100A
	Rabbit anti-Goat IgG-heavy and light chain Antibody Biotinylated	A50-100B
	Rabbit anti-Goat IgG-heavy and light chain Antibody FITC Conjugated	A50-100F
	Rabbit anti-Goat IgG-heavy and light chain Antibody HRP Conjugated	A50-100P
	Donkey anti-Goat IgG-heavy and light chain Antibody Affinity Purified	A50-101A
	Donkey anti-Goat IgG-heavy and light chain Antibody Alkaline Phosphatase Conjugate	A50-101AP
	Donkey anti-Goat IgG-heavy and light chain Antibody FITC Conjugated	A50-101F
	Donkey anti-Goat IgG-heavy and light chain Antibody HRP Conjugated	A50-101P

Catalog No.

A50-100AP

MATERIAL SAFETY DATA SHEET

Contents Description

This product contains Alkaline Phosphatase conjugated Rabbit immunoglobulin G (IgG) protein in 50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl₂, 0.1 mM ZnCl₂ containing 0.09% NaN₃.

Hazardous Ingredients

Rabbit IgG - No known toxicity for this biological material.

Bovine Serum Albumin (BSA) - No known toxicity for this biological material. (Applies only to buffers containing BSA.)

Alkaline Phosphatase - No known toxicity for this calf intestine material.

Sodium Azide (NaN₃) - AZIDE FORMS EXPLOSIVE CHEMICAL COMPOUNDS WITH LEAD AND COPPER PLUMBING. CARE MUST BE TAKEN TO WASH WASTE DOWN DRAINS WITH LARGE VOLUMES OF WATER.
LD50 oral mouse - 27 mg/kg.

Wash all affected areas with large volumes of water and if swallowed consult your physician immediately.

The above information is believed to be correct but does not purport to be all-inclusive and is intended to be used only as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state, or local laws or regulations.