

Biomedical Technologies Inc.

378 Page Street • Stoughton, MA 02072 USA • Phone: (781) 344-9942 Fax: (781) 341-1451 **Web:** www.btiinc.com

DATA SHEET

HUMAN FIBRONECTIN ANTISERUM

Catalog No: BT-551

Quantity: 1ml

**Protein
Concentration:** 8.14mg/ml

Lot No: 5511196

Packaging: BTI anti-Human Fibronectin IgG (rabbit) is packaged in 1ml aliquots with 0.1% sodium azide added as preservative.

Preparation: Rabbits were immunized with BTI Human Fibronectin (**BT-551**) which is approximately 95% pure (5% SDS acrylamide gel electrophoresis).

**Ochterlony
Analysis:** A single precipitin line obtained against human plasma (1 to 10 fold diluted) and a reaction of identity with purified fibronectin.

**Storage &
Storage:** This product contains sufficient NaN_3 for 4°C storage. Minimum one year.

FOR RESEARCH USE ONLY