

Biomedical Technologies Inc.

378 Page Street • Stoughton, MA 02072 USA • Phone: (781) 344-9942 Fax: (781) 341-1451 **Web:** www.btiinc.com

DATA SHEET

Rabbit anti-Human Fibronectin (Affinity Purified)

Catalog No: BT-551P

Quantity: 100ug

Lot No: 5511203

Concentration: 690ug/ml

Preparation: Obtained from immobilized human plasma fibronectin and dialyzed against phosphate buffered saline.

Ochterlony: Yields a single precipitin band against human plasma (1 to 10 fold dilutions) and reaction of identity with purified fibronectin.

Storage & Stability: Contains sodium azide. Storage at -20°C for 1 year.

FOR RESEARCH USE ONLY