



The Immunogold Specialist

Product Analysis

EM Bovine Serum Albumin: 20nm Gold

Product Information

Product Code: EM.BSA 20
Batch Number: ****
Protein Absorption: Not Applicable
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA
15mM NaN₃ (preservative) and 20% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	nothing detected	Passed
Optical Density:	4.1	Passed
Mean Diameter:	20.0nm	Passed
Coefficient of Variation:	3.7%	Passed
Percentage Single Particles:	98%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBinternational
Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK
Tel: 44 (0) 2920 747 232
Fax: 44 (0) 2920 747 242
E-mail: info@bbigold.com
www.bb-international.com

A complete catalogue of products can be found on our website or CD Rom