



The Immunogold Specialist

Product Analysis

EM Cationic Colloidal Gold: 5nm Gold

Product Information

Product Code: EM.CG5 5
Batch Number: ****
Protein Absorption: Not applicable
Protein Concentration: 18.1 µg / ml
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN₃ (preservative) and 20% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	<100pg Mucin	Passed
Optical Density:	4.2	Passed
Mean Diameter:	4.9nm	Passed
Coefficient of Variation:	5.4%	Passed
Percentage Single Particles:	Varies with changes of pH	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBinternational

Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: info@bbigold.com

www.bb-international.com

A complete catalogue of products can be found on our website or CD Rom