



The Immunogold Specialist

Product Analysis

LM Cationic Colloidal Gold: 5nm

Product Information

Product Code: LM.CGC 5
Batch Number: ****
Protein Absorption: Not Applicable
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA,
15mM NaN₃ (preservative) and 20% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected: 15pg after silver enhancement **Passed**
Optical Density: 3.1 **Passed**

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBinternational

Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: info@bbigold.com

www.bb-international.com

A complete catalogue of products can be found on our website or CD Rom