

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Canine Factor IX

Concentration: 0.75 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.133 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM
Benzamidine/ pH 7.4

Storage: < -60 °C

Expiration Date: 12/2006

Canine Factor IX is purified from canine plasma via ion and affinity chromatography. The protein is >90% as judged by SDS-PAGE.

Thank You,
Quality Assurance